

# Croyde Floor Mat

Space saving, robust & portable



- A rectangular shaped floor mat to support active movement at floor level during childbirth
- Innovative space saving design allows the mat to be folded in half for transport & storage
- Integrated grab handles are ergonomically positioned to allow the mat to be folded over and carried
- High friction base material is exceptionally durable and ensures positioning is stable



CROYDE MEDICAL

Maternity Equipment Specialists



- Lightweight and highly portable for movement around the birthing room
- Integrated handles are positioned to allow the carer to carry it without twisting
- Ease of movement minimizes the risks of inappropriate handling / dragging
- Designed to be wiped down in-between use in the same way as a hospital mattress
- Available to purchase with Croyde positioning pillows, birthing bags & wedges

## Specification

Length (Folded)	100cm
Length (Unfolded)	200cm
Width	120cm
Height (Folded)	12cm
Height (Unfolded)	6cm
Weight	4kg
Cover	Dartex Medical Grade
Crib	5
Chlorine Compatible	Yes, up to 10,000ppm
Warranty	1 Year



**CROYDE MEDICAL**  
Maternity Equipment Specialists

Croyde Medical, Unit 7, Regal Way,  
Faringdon, Oxon SN7 7BX