

Croyde Positioning Pillow

Comfortable, natural positioning support




- A triangular shaped beanbag delivering enhanced levels of comfort and support
- Offering two distinctive backrest angles
- 60° in the upright position
- 30° in the resting position
- Upright position supports back, neck, shoulders & head, with just one pillow
- Beads form around the body for optimum comfort & stability
- Resting position is ideal for examination / palpation / relaxing



CROYDE MEDICAL

Maternity Equipment Specialists



- Regular pillows are difficult to find, thin, easily compressed and are not steady enough to support sitting-up positions
- The Croyde Positioning Pillow delivers upright positioning with exceptional levels of comfort and stability
- Supports multiple, upright breast-feeding positions & reduces the core strain on the woman, especially after abdominal surgery
- Designed to be wiped down in-between use in the same way as a hospital mattress
- Available with faux suede, faux leather or cotton fabric covers for home use

Specification

Length	35cm
Width	50cm
Height	60cm
Weight	0.75kg
Cover	Dartex Medical Grade
Seams	Welded
Crib	5
Chlorine Compatible	Yes, up to 10,000ppm
Zip	Yes, covered with waterfall flap
Warranty	1 Year



CROYDE MEDICAL

Maternity Equipment Specialists

Croyde Medical, Unit 7, Regal Way,
Faringdon, Oxon SN7 7BX