

Croyde Birthing Bag

Comfortable, versatile & adaptable



- A barrel shaped birthing bean bag that helps the user find their most comfortable and supported positions during ante-natal, intrapartum & post-natal usage
- Helps support upright positioning and active movement during labour
- Can be shaped and moulded to suit the needs of the user
- Designed to be hugged, cuddled, & rested on



CROYDE MEDICAL

Maternity Equipment Specialists



- Immersive beads deliver exceptional levels of comfort with stability and support
- Promotes upright & active movement
- Can be used to support side lying positions with left lateral leg elevation
- Designed to be wiped down in-between use and is compatible with all chlorine-based solutions and wipes
- Zipped with a waterfall flap to prevent fluid ingress
- Available to purchase with Croyde floor mats, pillows & wedges

Specification

Height	65cm
Diameter	50cm
Weight	1.75kg
Top up beads	Optional
Cover	Dartex Medical Grade
Seams	Welded
Crib	5
Chlorine Compatible	Yes, up to 10,000ppm
Zip	Yes, covered with waterfall flap
Warranty	1 Year



CROYDE MEDICAL
Maternity Equipment Specialists

Croyde Medical, Unit 7, Regal Way,
Faringdon, Oxon SN7 7BX