

Oceanis Birthing Pool

Spacious, safe & personalized



- A beautiful birthing pool supporting freedom of movement and a range of options to enhance the birth experience
- Under-water, cool light LEDs enhance the examination process for the Midwife
- Birthing seat with perineal cut out and arm rests support upright resting positions
- Spacious base supports active movement
- Pop-up waste controlled from outside of the pool with a simple twist-release
- Can be installed up against an IPS wall or in a freestanding room set-up
- Exceptionally easy to clean with smooth surfaces that are easily accessible



CROYDE MEDICAL

Maternity Equipment Specialists



Additional Options to Personalize your Birthing Pool:

- External shell of the pool can be chosen in any colour to match / contrast the room design & deliver that spa-like feel
- Entry steps with handle lock into place using the unique flange system to make entry/exit significantly safer and easier
- Removable get-in aid can also support women entering the pool and provide an additional seat
- Grab handles can be fully integrated to the pool
- Underbody lighting can illuminate the floor area and be programmed to frequently change colour
- Bluetooth audio is available for women to connect to and play music through the shell of the pool



Specification

Length	190cm	Volume (to fill level)	450 litres
Width (Max)	132cm	Install Location	IPS or Free Standing
Height	74cm	Filling Taps / Mixers	Can be chosen to suit
Weight (Empty)	110kg	Integrated filling system	Available for free standing (DN19)
Weight (to Fill Level)	560kg	Warranty	2 Years



CROYDE MEDICAL
Maternity Equipment Specialists

Croyde Medical, Unit 7, Regal Way,
Faringdon, Oxon SN7 7BX