

# Croyde Positioning Wedge

Supporting women, midwives & birth partners



- A uniquely shaped foam wedge that offers a range of positioning options during childbirth
- The wedge is designed to support upright and active positioning during labour
- It can also be used as a comfortable seat by the woman, Midwife or birth partner
- The wedge is available as a stand-alone item but is supplied as standard with the Croyde Birthing Couch system



CROYDE MEDICAL

Maternity Equipment Specialists



- Lightweight and highly portable for movement around the birthing room
- Double covered in high quality Dartex fabric to reduce fluid ingress potential
- Visco-elastic foam top layer provides exceptional levels of comfort
- High resilience foam underneath provides stability and depth
- Designed to be wiped down in-between use the same as a hospital mattress
- Available to purchase with Croyde positioning pillows, birthing bags & wedges

## Specification

Length	60cm
Width	60cm
Height	43cm
Radius	50cm
Weight	4kg
Cover(s)	Dartex Medical Grade
Zip	Yes. Covered by waterfall flap
Crib	5
Chlorine Compatible	Yes, up to 10,000ppm
Warranty	1 year cover, 3 years foam



**CROYDE MEDICAL**  
Maternity Equipment Specialists

Croyde Medical, Unit 7, Regal Way,  
Faringdon, Oxon SN7 7BX