



CROYDE MEDICAL
Maternity Equipment Specialists

CROYDE GYNAECOLOGICALEXAMINATION COUCH

The Croyde Gynaecological Examination Couch is designed for thorough assessments whilst keeping the patient comfortable.



Designed for general gynaecology procedures with a high standard specification, available in a choice of colours.

The positioning options available create a comfortable environment for the woman, as well as an optimal working environment for the clinician.

The range of additional options includes a triple foot switch, a fold-down leg extension and a removable leg extension to help build the ideal procedural couch.



POSITIONING

The chair can be easily positioned with electronic adjustment of the backrest, hi-lo & trendelenburg tilt. The -15 degree tilt feature delivers perineal uplift for enhanced visual access.

Adjustable lithotomy poles enable safe and flexible positioning of the legs for women up to 240kg and the retractable wheelbase design delivers stability during usage.

The chair can also be laid flat.

STANDARD OPTIONS



STEEL FLUID BOWL



PAPER ROLL HOLDER



HEAD CUSHION

OPTIONAL EXTRAS



FOLD-DOWN LEG
EXTENSION



REMOVEABLE LEG
EXTENSION



TRIPLE FOOT
SWITCH

EVERYTHING YOU NEED TO KNOW

Click the links below to access the
Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Choices](#)

[CLICK HERE FOR INFORMATION ON HOW TO
BUY THE CROYDE GYNAECOLOGICAL
EXAMINATION COUCH](#)




CROYDE MEDICAL
Maternity Equipment Specialists

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX

