



**CROYDE MEDICAL**  
Maternity Equipment Specialists

## M-300 & M-400 BABY WEIGH SCALES

The M-300 and M-400 are a Class III MDD Approved portable baby scales, designed for accuracy, reliability and comfort.

The M-300 and M-400 Baby Weigh Scales are fantastically accurate and comfortable for the baby.

Both scales are super accurate down to the 5g mark. They also have Hold and Tare functions and adjustable feet, meaning they are perfect for any Birth Centres, Delivery Suites and Postnatal wards, as well as community usage.



### M-300

Weighing just 3kg, the M-300 the ultra-lightweight and portable version of the baby scale.

Accurate to just 2g up to 6kg with a maximum weight of 15kg.

Two AA batteries can power the M-300 for 600 hours.



### M-400

The M-400 is the bigger of the two options with a maximum weight of 20kg, despite only being 1kg heavier than the M-300. Comes with a digital height measurement rod.

Powered by six AA batteries with enough charge for 6,000 weigh-ins.



## OPTIONAL EXTRAS

The M-300 and M-400 can be fitted to the transport trolley to keep the scales safe and easy to move between rooms. The trolley comes with 2 brake wheels as standard for safety. Just order the scale with the trolley and we will mount them for you.

Both scales can also fit into a rucksack or carry case that are available as extras. There is also a power adapter that can power both scales.

### [CLICK HERE FOR INFORMATION ON HOW TO BUY THE M-300](#)

Click the links below to access the  
[Croyde Knowledge Hub](#)  
[Dimensions & Specifications](#)  
[Delivery & Implementation](#)  
[Cleaning Instructions](#)  
[Warranty & Servicing](#)

### [CLICK HERE FOR INFORMATION ON HOW TO BUY THE M-400](#)

Click the links below to access the  
[Croyde Knowledge Hub](#)  
[Dimensions & Specifications](#)  
[Delivery & Implementation](#)  
[Cleaning Instructions](#)  
[Warranty & Servicing](#)



**CROYDE MEDICAL**  
Maternity Equipment Specialists

## CONTACT US

Tel: 01367 241194  
Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)  
[www.croydemedical.co.uk](http://www.croydemedical.co.uk)

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX

