

# MATERNITY & NEONATAL EQUIPMENT DIRECTORY



  
**CROYDE MEDICAL**  
Maternity Equipment Specialists



Framework Ref: 2020/S 252-638069



Hi there,

I'm Stuart Hobin, Managing Director at Croyde Medical, and I'm thrilled to share this year's **Maternity and Neonatal Equipment Directory** with you.

This year's directory is our most comprehensive yet, featuring a wide range of new and innovative equipment.

In 2024, Croyde Medical proudly acquired **QED Scientific Ltd**, a UK-based manufacturer of warming cabinets and a distributor of maternity and neonatal equipment. As a result, we're excited to introduce several new products, including:

- QED blanket warming cabinet range
- QED fluid warming cabinet range
- Breast milk fridge range
- OKM-740 neonatal cot with integrated warming function

We've also added several other new and exciting products to our range:

- Croyde Baby Bath
- Croyde Twin Cot
- Aquatherm warming mattress & hot cot
- Vibwife birth mobility system (compatible with the AVE 2 birthing bed)
- Elios 3-in-1 mobile phototherapy and heating light

In addition, we've also been busy redesigning and improving some of the existing equipment.

As an example, we've changed the specification of the Croyde Mood Tube to include the new Croyde colour palette – calming hues that gradually transition over a 10-minute cycle. This feature helps empower the birthing person by allowing them to be more aware of contraction timings in sync with the colour changes.

We've also upgraded our website to show you more images and videos of the products you're keen to learn more about, to help you transform your maternity/neonatal environments with confidence.

As an NHS-approved supplier, all products in the directory are available through the NHS Supply Chain. And on top of that, we're proud to be a Carbon Neutral supplier!

We really hope you enjoy this year's directory. Thanks for being part of the Croyde journey 😊.

A handwritten signature in black ink, reading 'Stuart'.

Stuart

	Page Number	Maternity	Neonatal	Assessment	Theatres
Birthing Beds	2	✓			
Birth Mobility Systems	4	✓			
Feature Wallpaper	5	✓	✓	✓	
Birth Pools & Waterbirth	6	✓			
LED Sky Ceilings	10	✓	✓	✓	
Starry Sky Ceiling	10	✓			
Mood Lighting	10	✓	✓	✓	
Birthing Couch	12	✓			
Floor Mats, Beanbags & Pillows	12	✓			
Hooks, Ropes & Slings	14	✓			
Birth Stools & Balls	16	✓			
Low Working Kneeler	16	✓			
Co-sleeper Cot	18	✓			
Baby Changing	18	✓	✓		
Baby Assessment Table	19	✓	✓		
Height Adjustable Cots	20	✓	✓		
Cots and Neonatal Warming	22	✓	✓		
Twin Cot	24	✓	✓		
Phototherapy Light	24	✓	✓		
Baby Baths	25	✓	✓		
QED Warming Cabinets	26	✓	✓		✓
Breast Milk Fridges	27	✓	✓		
Maternity Recliners	28	✓	✓		
Recliner Chairs	30	✓	✓		
Fold-Away Beds	30	✓	✓		
Sleeper Chairs	31	✓	✓		
Wall Beds	32	✓	✓		
Double Beds	33	✓	✓		
Infant Feeding Chair	34	✓	✓		
Chairs & Overbed Tables	35	✓	✓	✓	
Examination Couches & Chairs	36			✓	
Mobile Examination Lights	37	✓		✓	



The AVE 2 brings the latest innovation to hospital delivery bed design. It has embedded patient controls, rotating side rails and folding grab handles. No extra storage or attachment of extra accessories is necessary.



MPC code: PPB-AX

## THE LATEST TECHNOLOGY

The AVE 2 birthing bed is now supplied with the upgraded Duocoat Premium mattress, the only birthing bed mattress with a visco-elastic memory foam top layer for enhanced comfort.

The Duocoast mattress has a double cover for protection against fluid ingress and is waterproof, breathable, zipped with welded seams and cleanable with hospital disinfectant solutions.





## SWIVELLING FOOT SECTION

The foot section of the AVE 2 never needs to be lifted off the bed. It simply lowers and stows under the bed in seconds, saving valuable time in emergency situations and avoiding the need to detach and store items.



## POSITIONING OPTIONS

The AVE 2 birthing bed delivers a plethora of birth positioning options to help women find their most comfortable profiles during childbirth. Upright and active postures are easy to achieve and are supported with comfortable positioning aids.

## CROYDE PEANUT BALL

Peanut balls have shown they can reduce labour time, regulate emotions and save your energy for when the birthing person is ready to push. This hypoallergenic BPA ball fits snugly between the legs to maintain optimal pelvic position.






## VIBWIFE

Vibwife is the only patented birth mobility system, replicating manual mobilisation techniques. Paired with the AVE 2 Birthing Bed, this mattress helps to support physiological & active birth using continuous oscillating movements. This movement of the mattress also helps with pain management and improves the birth experience.

MPC code: TBC







## FEATURE WALLPAPER

Bring colour to any hospital space using the Croyde Feature Hospital Wallpaper. A custom made wallpaper is an ideal way to cover a large area and make a great impact.

There is a wide range of designs available, from photographs of a local area to a calming coastal scene or landscape. These uplifting designs can transform clinical areas into welcoming spaces and help to create oxytocin-boosting birth environments.

Usually taking just a couple of hours to complete, this feature wallpaper can be installed without causing too much disruption to normal hospital routines and appointments.







## OCEANIS BIRTH POOL

MPC Code: Oceanis

The Oceanis birth pool is a warm, spacious pool delivering freedom of movement in water. It is supplied as standard with LED underwater lights, externally controlled pop-up waste and integrated birth seat.

The innovative “lock-in” entrance steps with handrail make entry and exit significantly easier and safer for the woman. They lock to the pool to give guaranteed stability when in use and move away on integrated wheels when not.

The external colour of the pool can also be specified to suit any room décor and under-pool lighting is available to introduce soft mood lighting to the room.





## COMBIPOOL

The CombiPool is a luxury birth pool that delivers an exceptional water birth experience. The shape and aesthetics are soft and gentle, contributing to a calming birthing environment.

Non-slip steps and an angled handrail provide excellent stability when entering and exiting the pool. The internal floor space is large enough to support freedom of movement and offers several pre-defined positions.

The design of the sides achieves the optimum working height for the attending professionals, minimising manual handling risks. The shaped end also allows the midwife or birth partner to get closer to the water by placing their legs underneath the rim.



MPC Code: COMBI-Pool

## GENEROSA

The Generosa Birth Pool is a compact yet spacious birth pool, making it well-suited to smaller birth rooms that are tight on space. Designed with an ergonomic built-in seat as well as integrated handles and foot supports for comfort and support, the Generosa Birth Pool creates an optimal environment for the birthing person.

MPC code: CAP-GBP01



## PELVITUB

The Pelvitub is a fully submersible mat that can be used for kneeling, sitting or squatting in birthing pools. It is completely waterproof and easy to wipe down after use. Suction cups ensure that the mat doesn't float and it is lightweight and comfortable to use.

MPC code: PELVIK



# CROYDE UNIVERSAL BIRTH POOL STEPS

Croyde Universal Birth Pool Steps have been thoughtfully created with both the birthing person and hospital staff in mind, creating a safe and fit for purpose solution. They can be safely used with any birth pool.



## SAFER BIRTH POOL ENTRY/EXIT

Designed to be used safely with any birth pool, the double step design provides easy access to the pool and the handrail gives additional support & reassurance.



## STABILITY & BALANCE


The Croyde Stability Bar has been added underneath the base to ensure stability. This extra weight lowers the centre of gravity, making it secure and tip-resistant.

The base is also flared to maximise the surface area.



## MOVEMENT & STORAGE

The folding handle means that the steps can be folded to the size of a small suitcase and easily stored. This design means the steps can be either pushed or pulled into position. They can also be tipped back onto their wheels to be moved around.



## VERSATILE AND CLOSE

The uniquely shaped handrail allows the steps to be placed closer to the pool, reducing the amount of stretching required to reach the rail.

Thoughtfully chosen, the smaller circumference of the rail means that it is comfortable to hold, even for small or swollen hands.

## CROYDE LED SKY CEILINGS


Bringing the outside – inside. LED Sky Ceilings help to create a more expansive, relaxed environment. Designed to help reduce anxiety and promote feelings of calmness and warmth, they can be installed into any tiled ceiling.

They are popular in pool rooms, bereavement suites, waiting areas and rooms with little or no natural light.

Available in a wide range of images and sizes.



MPC code: LED SKY



MPC code: ROUND / AVOCADO

## CROYDE STARRY SKY CEILINGS

A range of ceiling shapes and designs to enhance the room appearance. Twinkling stars and LED edge lighting create a whole new environment when the lights are low.

Excellent for use in pool rooms and can be tailor-made to match the exact shape of the pool.



## CROYDE TEN MINUTE MOOD TUBE

A mood enhancing column of light that cycles through a curated colour palette on a ten-minute cycle.

With each colour change, the birthing person can feel empowered to recognise their contraction timings, and be provided the opportunity to gain a better understanding of their labour progress.

Used as a focal point in the room, it can help to create a more relaxed and soothing environment.

MPC code: CHROMA


## CROYDE LIGHT CRYSTAL

The Croyde Light Crystal cycles through a magical world of colours to deliver a calming, sensory experience.

It is placed in the corner of the room and patterns of light are projected through the crystal dome to produce a mesmerising, light display.

The crystal works straight out of the box and just needs to be plugged in with no ongoing maintenance.

Perfect for areas striving to deliver relaxed environments.



MPC code: DOME



## CROYDE AIR TUBE

Experience the captivating display of swirling beads as they spiral through the Croyde Air Tube.

The beads in the tube are propelled upwards and cascade down, creating a beautiful vortex pattern.

A similar effect to bubbles in a bubble tube but without the necessary maintenance!

MPC code: PASSIVE



## CROYDE® BIRTHING COUCH

The innovative Croyde Birthing Couch is designed to support active and natural childbirth, helping women to find their most comfortable positions during labour and delivery. The Croyde Birthing Couch can be used as a full-length couch, folded over and used to support standing labour or used as a seat; and even be used as a double bed for post-natal resting.

Double covered in high quality, medical grade, breathable fabric, it is resistant to fluid ingress, easy to clean and compatible with hospital wipes and chlorine based cleaning solutions.

MPC code: CBC



## CROYDE® POSITIONING PILLOW

A triangular shaped beanbag that comfortably supports women when sitting upright or resting. Delivering a 60° backrest in the upright position and 30° in the lower position, it shapes and moulds to the body to enhance comfort whilst avoiding having to find lots of regular pillows.

Compatible with hospital wipes and chlorine based cleaning solutions.



MPC code: CPP

## CROYDE® BIRTHING BEANBAG

A barrel shaped beanbag designed to retain its shape and support upright positions. The bag can be shaped and moulded to the specific needs of the woman and is comfortable to cuddle and rest upon.

Compatible with hospital wipes and chlorine based cleaning solutions.



MPC code: CBB

# CROYDE® POSITIONING PRODUCTS

When used together, the Croyde products offer a plethora of positioning options to help the woman, Midwife and birth partners to become engaged.



## CROYDE® FLOOR MAT

A large, comfortable floor mat delivering space for active movement during labour. It can be folded over to be stored in a compact, half-length position and is easily carried using the integrated handles.

Compatible with hospital wipes and chlorine based cleaning solutions.

MPC code: CFM



# MULTITRAC BIRTH SYSTEM


The Multitrac Birth System promotes close and active labours between the birthing person and partner. It is supplied with a plug-in Dullstein Birth Stool, Pelviball with fixed base, foot mat for comfortable kneeling and integrated rope system for vertical hanging. Multitrac is a complete system for supporting upright & active birth.

**“Since purchasing the Multitrac earlier this year, women have been able to adopt positions that support physiological birth. The Multitrac enables the birth partner to provide support and comfort during labour and birth. It is easy to use, and the integrated rope sling has been a great tool for women.”**

Abigail Holmes, Deputy Head of Midwifery, Cardiff and Vale University Health Board



MPC code: MTRAC




## ROPE SLINGS

Are supplied with the Multitrac as standard but are also available to purchase individually in a range of colours.

MPC code: SDIDY







## VARIOTRAC® SWIVEL SYSTEM

The Variotrac® swivel system is an integrated system that allows women to use an overhead rope in whichever position they like.

It has 360° swivel capability and locks at 15° increments so you can position it above a birthing pool, couch, bed, birth stool, chair or any other platform on demand. Rope slings are safely and easily attached onto any of the 5 hook positions and are available in a range of colours.

Variotrac can be mounted into the ceiling, wall or floor to suit the room.

MPC code: VARIO D (Ceiling)

MPC code: VARIO S (Floor)

MPC code: VARIO W (Wall)

## STATIC CEILING HOOK

Our ceiling-mounted hooks can be installed into any birth room ceiling to give women the option to vertically hang, using a rope sling and adopt their most comfortable positions.

MPC code: STAT-HAK



## DULLSTEIN® BIRTH STOOL

The Dullstein® Birth Stool is designed to support upright birth and is covered with a soft imitation leather cover to enhance comfort during usage. It has sealed seams to prevent fluid ingress and is available in a choice of colours – yellow, blue and orange with height options of 33cm or 38cm.

The Dullstein® can also be combined with the perfectly shaped inlay mat to support additional active birth positions. It is lightweight, easy to wipe clean and has a max SWL of 135kg.

MPC code: DGH - DGH-alu



## CROYDE KAYA BIRTH STOOL

Strong, stable and easy to clean, this innovative birth stool is able to support a range of positions. It can be used under water for additional options in birthing pools and can also be used for rocking.



MPC code: KBS







MPC code: ERG001

## ERGOKNEELER

The innovative Ergokneeler is an ergonomic kneeling system that protects knees and promotes good musculoskeletal health.

It can enhance comfort levels for Midwives who are supporting women during active labour and helps to reduce the compression forces and postural discomfort associated with low working height positions.

Moulded in polyurethane with active antibacterial properties, the ergonomic design offers supported kneeling and sitting positions whilst also being lightweight, portable and easy to clean.



## PELVIBALL

The Pelviball can be used on its own or in conjunction with the Multitrac or Duoetto birthing systems.

The ball is attached to a base plate to offer maximum stability and security during labour, whilst also preventing it from inadvertent rolling. The ball is easy to wipe clean, is anti-burst and easy to top up with air over time.

Optional covers are available with disposable and washable versions.

MPC code: P-ball(65) / P-ball(75)

MPC code:  
BBB(65) / BBB(75)

## CROYDE BIRTHING BALLS

Croyde birthing balls feature anti-burst technology and material that is 5 times thicker and more durable than standard gym balls.



## CROYDE CO-SLEEPER COT

The Croyde Co-Sleeper Cot is shaped to give immediate access to the baby with the safety of lipped edges and integrated side rails. It is designed for use in hospitals and helps to deliver a home-from-home environment.

MPC code: COSLC

## CROYDE BABY CHANGE STATION

The Croyde Baby Change Station has clinical application with a homely appearance. It's ideal for Birth Centres and Midwife-Led Units. The changing unit features a smooth tabletop surface with a curved upstand, deep, lockable drawers for safe storage, a soft oak finish and integrated 75mm castors with brakes for ease of movement.

The soft-close drawers provide enough space to store essential items. For example, you can store basic neonatal resuscitation equipment to replicate the home birth setup.

Standardised equipment is then accessible, but not visible, a scenario that is well-received during diligence reviews.

MPC code: CBS-M







## HEIGHT ADJUSTABLE BABY ASSESSMENT TABLES

Available in a standard or compact size, the Croyde Height Adjustable Baby Assessment Table allows clinicians to position the table at the perfect height for neonatal examinations and baby changing.

Depressing the foot pedal enables staff or parents to adjust the table to the ideal position and protect them from physical stress and strain. NIPE screening becomes instantly more comfortable and stress-free.

The high-height position makes cleaning easier, while the low-height position makes the table accessible for all heights of carers, including wheelchair users. It is lightweight and easy to move between bed spaces, with brakes to secure its location.

Choose from a storage cupboard with shelving or a drawer unit. Both units have soft-close doors and ample space for all items.



Standard MPC code: SPTM-02-D-HA  
Compact MPC code: SPTM-02COMPACT-D-HA



## CROYDE HEIGHT ADJUSTABLE OVERBED COT

The Croyde Height Adjustable Overbed Cot glides under the bed, putting the baby in front of Mum. It allows constant and easy access, facilitating more frequent skin-to-skin contact, and promotes maternal bonding.

The mother can place the baby back in the cot without leaving the bed, reducing the risk of her falling asleep with the baby because the cot is too far away.

The Overbed Cot can be raised and lowered in one simple smooth movement to make access to the baby easier for women who have undergone caesarean sections or complicated births.

The open design makes it easy to clean, and a wide range of colours are available to suit any room theme.

MPC code: BM/04/GS-KR-01ST







## CROYDE HEIGHT ADJUSTABLE COT + STORAGE

Chosen by birth centres, postnatal areas, NICUs, SCBUs and HDUs, the Croyde Height Adjustable + Storage cot offers the convenience of customisable and practical under-cot storage.

The large cupboard is accessible from both sides of the cot with double soft-close doors, ensuring it is suitable for small rooms or tight spaces.

Height adjustment enables the cot to be raised or lowered to the desired height to limit bending or stooping, reducing musculoskeletal strain on the midwife's back and shoulders.

You can specify an optional foldable shelf or oxygen cylinder holder, and a wide range of colours are available to suit any room theme.

MPC code: BM/05/GS



## AQUATHERM WARMING MATTRESS/HOT COT

The Aquatherm warming mattress provides heat and warmth to regulate a neonatal patient's temperature. It uses 'Aqualgel' technology to deliver a non-water based warming system.

It incorporates a heating element for safe, evenly spread heat and two highly-sensitive sensors for fail-safe temperature control. The control unit has user-friendly settings to regulate the heat output of the mattress and also contains an alarm system for enhanced patient safety.

The Aquatherm is available as a stand-alone mattress and control unit, or as a hot cot solution with the Croyde Height Adjustable Cot. When used as a hot cot system, the swivel bracket provides immediate access to the control unit and discreetly locates it inside the frame to avoid accidental collisions.

Access holes are created in the crib's basin to allow the cable to conveniently connect to the mattress without kinking. The shape also fits perfectly within the crib to minimise gaps. When disconnected for transport, Aqualgel technology ensures that residual heat is retained for longer.



MPC code: Aquatherm





## CROYDE OKM-740 BABY WARMING BED

An electrically adjustable neonatal cot with integrated warming mattress and LCD display controls, suitable for single or twin occupancy.

The OKM-740 has an integrated heating system and gel mattress with temperature settings between 30°C & 38°C in 0.1°C increments.

Skin and mattress probes accurately present temperature on the LCD colour display which also shows the SpO<sub>2</sub> reading and trend information.

The sides of the cot create a draught-free environment whilst still allowing ventilation. Each side can be dropped to allow access to the baby / babies as required.

The optional lid can be added to maintain optimal warmth, and can also be folded back to give access to the baby.

The cot also has a Trendelenburg function which can be selected at the push of a button on the control unit.

## TWIN COT

The Croyde Twin Cot features a spring-assisted height adjustment pedal, allowing for effortless raising and lowering in one smooth motion. This allows clinicians to adjust the height for neonatal examinations and mothers easy access to their newborns.




MPC code: SBS-2.0-Twin-12



The easy-to-use release buttons on each side allow the side panels to be lowered for convenient access.

Designed with safety in mind, both catches must be released before lowering each side panel.

When lowered, the side panels provide full access along the side of the cot and enable it to be positioned next to the bed, allowing mothers to safely reach their babies.



## ELIOS PHOTOTHERAPY LIGHT

ELIOS seamlessly integrates intensive phototherapy with heat treatment, providing a comprehensive solution for neonatal care. The lamp is designed for safe and easy placement over cribs and changing tables.

This versatility is facilitated by its robust, height-adjustable steel stand and an anti-tipping base equipped with swivelling wheels for smooth manoeuvrability.

MPC code: 8037



# CROYDE BABY BATH

The Croyde Baby Bath is a height adjustable baby bath with our Easy-Pour Tub for quick and simple emptying.

The tub can be raised and lowered using the foot pedal to tailor the height to suit any care giver and improve ergonomics when using or emptying.



The Easy Pour Tub design features a pointed spout to make emptying the bath efficient whilst reducing spillages. The handle is easy to grip when tipping and can also be used as a small shelf to store soap sachets, cotton wool etc when in use.

MPC code: CAP-TBB01

# TIPPER TUB

This innovative mobile baby bath enhances caregiver ergonomics and eliminates the need to manually lift and carry heavy bowls of water.

The integrated lifting column allows the tub to be raised to the height of the sink or drain where it is then tilted without any hand contact to empty without strain. The bowl is ergonomically shaped, heat-insulated to retain the warmth of the water and one-piece to make cleaning exceptionally easy.

Available in a choice of colours.



MPC code: TT-24000

## COUNTER CABINET: BLANKET WARMING

The QED Blanket Warming Cabinets are available in a range of sizes accommodating 4-50 cotton blankets in baskets. These cabinets can be used in range of settings, including birth environments, postnatal & neonatal, operating theatre recovery, critical care, A&E and other hospital wards where rapid patient warmth is important.

All QED Warming Cabinets are equipped with PID digital temperature controls. Blanket warmers have a maximum temperature of 80°C but can be set to a desired temperature.

The chamber is made from medical grade stainless steel and the external of the cabinet is powder painted to allow for easy cleaning.



MPC code: BW-ECO



MPC code: FW-ECO

## COUNTER CABINET: FLUID WARMING

The QED Fluid Warmers are available in a range of sizes accommodating 9 to 240 one litre bottles of fluid.

Smaller cabinets are designed to fit underneath or on top of work surfaces. Larger cabinets are equipped with wheels for ease of movement.

Fluid warmers have a maximum temperature of 60°C but can be set to a desired temperature.

Each cabinet is fitted with a stainless steel interior for easy cleaning.

The stainless steel interior and powder-coated exterior are compatible with hospital-approved wipes and chlorine-based cleaning solutions.

## BREAST MILK FRIDGES

The Croyde Breast Milk Fridges are available in a range of sizes, from sixty litre capacity to four hundred litre capacity.

The digital control panel allows for the temperature selection, between 1-5 °C, with a single press 3 °C, setting for the optimum breast milk temperature.

Audible and visual temperature alarms immediately alert healthcare practitioners of any major temperature fluctuations and help to ensure the quality of the milk.

MPC code: BWR125







## COMBI CABINETS

The QED Combi Warmers come in a wide range of sizes dictated by the individual size of the fluid warmer and the blanket warmer.

These sizes do not have to be equal, meaning your chosen Combi Warmer can have more space dedicated to the fluid warmer or blanket warmer.

Fluid warmers compatible with the combi warmers range hold 15 to 72 one litre bottles of fluid.

Blanket warmers compatible with the combi warmers range hold 4 to 24 blankets.

Product code: FW/BW

## FLOOR STANDING CABINETS: FLUID AND BLANKET WARMING



Product code: FW ECO



Product code: BW ECO

# CROYDE MATERNITY & NEONATAL RECLINER

A multi-use bedside package, perfect for maternity wards. The recliner facilitates ideal infant feeding positions for mother and baby. The reclining backrest allows babies to feed in the semi-recumbent chest-to-chest position that is most comfortable for the mother. Padded armrests give comfortable arm placement and support the chosen feeding position. For partners and carers, the recliner positions them perfectly for skin-to-skin contact when bonding with the baby.

MPC code: ASTCMR



The padded foot section of the Croyde Maternity Recliner gently raises the level of the feet closer to the heart, improving lymphatic drainage, and reducing swelling in the feet, ankles, and legs after birth. By elevating the feet, the lower back also benefits from reduced strain.

The chair also allows overnight stays for partners without needing additional equipment, replicating the home environment and improving bonding with birth partners.

The Croyde Infant Feeding Pillow and Footstool are also available to purchase separately.

Available in a range of colours.





# CROYDE ELECTRIC MATERNITY RECLINER

The Croyde Electric Maternity Recliner Chair aids sitting, reclining, and rising at the touch of a button, and supports women who have undergone invasive surgeries. It makes the postnatal period more comfortable for women who have restricted mobility or have undergone a caesarean section by reducing the use of abdominal muscles and back strain by gently lifting them to a standing position.

For partners and carers, the chair positions them perfectly for skin-to-skin contact when bonding with the baby.

The Croyde Electric Maternity Recliner Chair provides a multi-use bedside package for maternity wards. It supports overnight stays without needing additional equipment, replicating the home environment and improving bonding with birth partners. The high-winged backrest gives neck and shoulder support and provides privacy in bay areas.

MPC code: CEMR-20/CEMR-25



# CROYDE RECLINER CHAIR

The Croyde Recliner Chair delivers fantastic versatility for patients and carers without taking up additional floor space - its footprint is similar to a regular bedside chair.

The chair provides fantastic comfort and support and effortlessly reclines to a lying position with a single, smooth action. This flexibility allows visitors to sit or rest at a patient's bedside without needing additional furniture, even during overnight stays.

The highly durable material is designed for high-traffic environments and it is compatible with hospital-approved wipes and chlorine-based disinfectant solutions.

All our hospital recliners are lightweight with larger castors than other chairs, making them easy to move around for cleaning and reorganising. The high backrest also means they can be moved around without bending over. The brake option ensures they can be safely locked when stationary.

Available in a range of colours.

# CROYDE FOLD-AWAY BED

The Croyde Fold Away Bed is a full-length single Bed with a comfortable mattress, easily deployed in two movements. It is durable and designed to be used on high traffic wards.

Lightweight and easy to move, it stores easily in the compact, folded position.

The solid bed deck means that cleaning the bed is simple, with no mesh bases, springs or difficult to clean areas. The mattress is also waterproof and can be cleaned using hospital approved cleaning solutions and wipes.



## LIGHTWEIGHT AND PORTABLE

The Croyde Fold Away Bed is easily moved into position by one person.



## TWO-STEP SET-UP

Fold down one side, then fold down the other. No spring-loaded mechanisms!



## TWO-STEP FOLD AWAY

Fold up one side, then fold up the other. Lockable castors keep the bed stable.





MPC code: STARGAZER



## CROYDE SLEEPER CHAIR

Perfect for use in hospital settings where overnight stays for partners/dependants are required. The Croyde Sleeper Chair is easily deployed from chair to comfortable, full-length bed in just a few seconds. It also features a chaise longue position. Integrated wheels make moving simple and lightweight and it is designed to be wiped-down between uses in the hospital setting.

Available in a range of colours.

**“The sleeper chairs are both spacious and comfortable, whilst not taking up too much room in our side rooms. They have helped us to ensure that our patient’s relatives have a better experience during overnight stays, which are often in very difficult circumstances”**

Inpatient Matron at Basildon & Thurrock University Hospitals NHS Foundation Trust

# CROYDE WALL BED

The Croyde Wall Bed is available as a single bed or a double bed and can be easily raised and lowered with just one hand.

The beds can be purchased independently or with alongside cabinetry such as wardrobes, cupboards, shelves or storage boxes.


Integrated lighting, USB charging points and power options are optional and a comprehensive delivery and installation service is available. The cabinetry can be adapted to accommodate gas systems, nurse call systems or other fixed wall systems.

A wide range of colours and finishes is available and each bed is supplied with a hospital-specification mattress that fits perfectly to the frame and easily folds away with the bed.


All items can be cleaned with hospital-approved wipes and chlorine-based solutions.




## SINGLE BED



## DOUBLE BED



## DOUBLE WITH CABINETRY







Above and above right: Bespoke installation of the Croyde Wall Bed, Croyde LED Sky Ceilings and Feature Wallpaper at The John Radcliffe Hospital, Oxford



MPC codes:  
Double wall bed: CDWB  
Double wall bed & cabinetry: CDWC+C  
Single wall bed: CSWB



## CROYDE BEREAVEMENT SUITE/ELECTRIC DOUBLE BED

Our electric profiling bed features hi-lo, backrest and thigh section adjustment via a back-lit handset control. It is available in 120cm or 140cm widths, or as a 90cm single bed if space is tight.

The bed has a weight capacity of 315kg and has integrated CPR release for emergency situations.

It is supplied with a hospital-specification mattress and can be cleaned with approved wipes and disinfectants. Choose from a range of lacquered wood finishes to create a home-from-home environment.



MPC code: XB1P8ABD106

## CROYDE STATIC DOUBLE BED

Our static double bed is supplied with a high headboard and low footboard to give the bed a homely aesthetic.

The lacquered wood finish enables the bed to be cleaned with hospital-approved wipes and disinfectant solutions.

It is also supplied with a hospital-specification mattress.

A range of colours and finishes is available.



MPC code: ABIBED-STDBL

MPC code: 365900



## CROYDE MIMMAMÀ INFANT FEEDING CHAIR

Specifically designed for infant feeding to accommodate a variety of different positions including cradle, cross and rugby ball options.

Ideal for twin feeding due to the dual backrest spaces. There is excellent support for the lower back, arms and shoulders and available with the optional footrest.

Available in a range of colours.



MPC code: E600



## MATERNITY CHAIR

Croyde maternity ward chairs are a fantastic addition to any unit and are ideal for use in post-natal wards and birth rooms alike.

They have exceptionally comfortable visco-elastic foam seating to enhance the user's sitting/resting experience and a stylish, supportive design.

Available in a wide range of colours, they can be ordered to suit any room décor. The vinyl fabric is also infection control compliant and designed to be wiped down in-between use.

Manufactured in the UK to exacting quality standards, the frames are stable, robust and provided with a 5 year warranty.

## VISITOR CHAIR

Croyde visitor chairs are comfortable, practical and can be colour coordinated to patient chairs and lockers.

They are available with a choice of fabric colours to complement modern room designs and colour schemes; whilst also having comfortable foam seating to enhance the user experience.

Manufactured in the UK using only the highest quality materials, they are strong, durable and designed for use in the hospital setting.



MPC code: PT10K/PT11K



MPC code: F833

## OVERBED TABLE

Croyde Overbed Tables are height adjustable and designed to allow the table to be used equally well with a bed or chair with its U-shaped base.

They are lightweight, easy to move on wheels, durable for longevity of use and easy to clean.

Table tops have lipped edges to prevent items sliding off and are available in blue, grey or beech finishes.

Tilting or non-tilting top options are available.

## CROYDE EXAMINATION COUCHES

A range of electric and hydraulically adjustable examination couches, all with excellent low height position for patient entry/exit.

Choose from 2 section or 3 section models and options are available for wheels or feet to suit the purpose.

Available in a range of colours.



MPC codes: SM2560 / SM3560

## CROYDE SUTURING COUCH

The Croyde Suturing Couch features electric positioning for height, backrest angle and tilt position.

It is supplied as standard with adjustable lithotomy leg rests, matching head support cushion, side support loops, debris tray and paper roll holder. An optional leg extension plate is also available.

The 70cm platform width enhances comfort during usage and the large castors and centrally operated brakes make it lightweight and easy to transfer between rooms.



MPC code: NV8683



## CROYDE OPERATOR AND SADDLE STOOLS

Height adjustable operator stools and saddle stools. Comfortable seating with options for wheels, feet and brake-on seating.

Available in a range of colours.



MPC codes: MC6164 & MC6144

## CROYDE DIASYS EXAMINATION LIGHT

A superb quality LED examination light capable of delivering high-level illumination for minor procedures.

Hands free brightness adjustment and automatic switch-on when you pull the light head into position, it is easy and practical to use when suturing.

Turns off after 30 minutes inactivity to avoid unnecessary usage.



MPC code: DIASYS-MOB-S0-UK

## CROYDE MINOR OPS E1 LIGHT

This examination lamp features high-quality, energy-efficient LEDs and a flexible extension arm to deliver illumination for minor procedures.

The mobile base of the light is stable, well balanced and durable with a radius footprint of just 350mm to make it easy to store.

The light head is available with a mobile base or can be mounted to a wall or ceiling.



MPC code:  
MAR-LED E1

SCAN THIS CODE BELOW TO VISIT THE  
EQUIPMENT WEBSITE AND VIEW EACH  
ITEM IN MORE DETAIL:



ALTERNATIVELY, YOU CAN VISIT:  
[WWW.CROYDEMEDICAL.CO.UK](http://WWW.CROYDEMEDICAL.CO.UK)

If you would like to discuss the specifications or suitability of  
any of the equipment in this directory, please call our helpline  
on 01367 241194 for further assistance.



CROYDE MEDICAL

Maternity Equipment Specialists

## CONTACT US

Tel: 01367 241194  
International Tel: 00 44 1367 241194  
Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)  
[www.croydemedical.co.uk](http://www.croydemedical.co.uk)

Purchase orders should be sent to: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)  
or faxed to: 01367 241657

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon SN7 7BX

