

MATERNITY EQUIPMENT DIRECTORY

All your maternity equipment needs, expertly supported.



CROYDE MEDICAL



CERTIFIED BY AFOR UK LIMITED



Medical Healthcare
Furniture
In partnership with
NHS Supply Chain

AVE 2 BIRTHING BED

Product Code: PPB-AX

NO LIFTING OR REASSEMBLY



The AVE 2 brings the latest innovation to hospital delivery bed design. It has embedded patient controls, rotating side rails and folding grab handles. No extra storage or attachment of extra accessories is necessary.

THE LATEST TECHNOLOGY

The AVE 2 birthing bed is now supplied with the upgraded Duocoat Premium mattress, the only birthing bed mattress with a visco-elastic memory foam top layer for enhanced comfort.

The Duocoat mattress has a double cover for protection against fluid ingress and is waterproof, breathable, zipped with welded seams and cleanable with hospital disinfectant solutions.


POSITIONING OPTIONS

The AVE 2 birthing bed delivers a plethora of birth positioning options to help women find their most comfortable profiles during childbirth. Upright and active postures are easy to achieve and are supported with comfortable positioning aids.



SWIVELLING FOOT SECTION

The foot section of the AVE 2 never needs to be lifted off the bed. It simply lowers and stows under the bed in seconds, saving valuable time in emergency situations and avoiding the need to detach and store items.



CROYDE PEANUT BALL

Product Code: BCS-M

Peanut balls have shown they can reduce labour time, regulate emotions and save your energy for when the birthing person is ready to push. This hypoallergenic BPA ball fits snugly between the legs to maintain optimal pelvic position.



VIBWIFE BIRTH MOBILITY SYSTEM

Product Code: D9U001VBB-0101



C-section rates are rising. Interventions are increasing. Epidurals are becoming the norm, whilst positive birth stories are becoming rarer. Evidence shows that active labour improves delivery outcomes, but when large portions of labour are spent on the bed this is almost impossible to achieve, especially when epidurals are used.

The Vibwife Birth Mobility System changes this. An innovative active mattress, inspired by the principles of Spinning Babies, that delivers in-bed mobility and replicates natural movement. Improving birth experiences and enabling active births for all.

Rocking Baby

A gentle motion to create asymmetry in the pelvis.

Sleeping Beauty

Changing the angle of the pelvis through a simple side to side motion

Infinite Harmony


Improving pelvic flexibility through a figure eight motion similar to walking or stair climbing.

Vibration & Jiggle

Producing calming effect to relax muscles for pain relief.



THE ACTIVE MATTRESS FOR BETTER BIRTHS



TECHNICAL SPECIFICATIONS

Specification	Value
Safe Working Load	220 kg
Width	855 mm
Weight	26 kg



OCEANIS BIRTH POOL

Product Code: Oceanis

The Oceanis birth pool is a warm, spacious pool delivering freedom of movement in water. It is supplied as standard with LED underwater lights, externally controlled popup waste and integrated birth seat.

The external colour of the pool can also be specified to suit any room décor and under-pool lighting is available to introduce soft mood lighting to the room.


The innovative “lock-in” entrance steps with handrail make entry and exit significantly easier and safer for the woman. They lock to the pool to give guaranteed stability when in use and move away on integrated wheels when not.



LED
COLOUR
CHANGE



LOCK IN
STEPS



COMBIPOOL BIRTH POOL

Product Code: COMBI-Pool

The Combipool is a luxury birth pool that delivers an exceptional water birth experience. The shape and aesthetics are soft and gentle, contributing to a calming birthing environment.

Non-slip steps and an angled handrail provide excellent stability when entering and exiting the pool. The internal floor space is large enough to support freedom of movement and offers several pre-defined positions.

The design of the sides achieves the optimum working height for the attending professionals, minimising manual handling risks. The shaped end also allows the midwife or birth partner to get closer to the water by placing their legs underneath the rim.



GENEROSA BIRTH POOL

Product Code: CAP-GBP01



The Generosa Birth Pool is a compact yet spacious birth pool, making it well-suited to smaller birth rooms that are tight on space. Designed with an ergonomic built in seat as well as integrated handles and foot supports for comfort and support, the Generosa Birth Pool creates an optimal environment for the birthing person.



CROYDE UNIVERSAL BIRTH POOL STEPS

Product Code: CUBPS



Croyde Universal Birth Pool Steps have been thoughtfully created with both the birthing person and hospital staff in mind, creating a safe and fit for purpose solution. They can be safely used with any birth pool.



VERSATILE AND CLOSE
 The uniquely shaped handrail allows the steps to be placed closer to the pool, reducing the amount of stretching required to reach the rail.
 Thoughtfully chosen, the smaller circumference of the rail means that it is comfortable to hold, even for small or swollen hands.



SAFER BIRTH POOL ENTRY/EXIT
 Designed to be used safely with any birth pool, the double step design provides easy access to the pool and the handrail gives additional support & reassurance.



STABILITY & BALANCE
 The Croyde Stability Bar has been added underneath the base to ensure stability. This extra weight lowers the centre of gravity, making it secure and tip-resistant. The base is also flared to maximise the surface area.



MOVEMENT & STORAGE
 The folding handle means that the steps can be folded to the size of a small suitcase and easily stored. This design means the steps can be either pushed or pulled into position. They can also be tipped back onto their wheels to be moved around.

PELVITUB

Product Code: PELVIK

The Pelvitub is a fully submersible mat that can be used for kneeling, sitting or squatting in birthing pools. It is completely waterproof and easy to wipe down after use. Suction cups ensure that the mat doesn't float and it is lightweight and comfortable to use.





LED SKY CEILINGS

Product Code: LED SKY

Bringing the outside – inside. LED Sky Ceilings help to create a more expansive, relaxed environment. Designed to help reduce anxiety and promote feelings of calmness and warmth, they can be installed into any tiled ceiling.

They are popular in pool rooms, bereavement suites, waiting areas and rooms with little or no natural light.

Available in a wide range of images and sizes.

STARRY SKY CEILINGS

Product Code: ROUND/AVOCADO

Croyde Starry Sky Ceilings are designed to help midwives transform maternity environments with bespoke low-level lighting.

Ceilings are available in a range of shapes and designs to enhance the room appearance.

The midnight blue back panel, twinkling stars and LED edge lighting create a whole new environment when the lights are low.

Excellent for use in birth rooms and birth pools, they create a calm, relaxing space. They can also be tailor-made to match exact pool shapes.



CROYDE TEN MINUTE MOOD TUBE (PLUS AUDIO)

Product Code: Chroma

A mood enhancing column of light that cycles through a curated colour palette on a ten-minute cycle.

With each colour change, the birthing person can feel empowered to recognise their contraction timings, and be provided the opportunity to gain a better understanding of their labour progress.

The Croyde Ten Minute Mood Tube is also available with integrated Bluetooth speakers for easy connection to a mobile device, allowing the birthing person to play a calming audio of their choice for a more soothing experience.

Ten Minute Mood Tube plus Audio - Product Code: CT100P15-MEDx3-BA



TEN MINUTE COLOUR CYCLE



OPTIONAL INTEGRATED BLUETOOTH SPEAKERS

Unlike standard off-the-shelf sensory light solutions that are designed to be stimulating, the Croyde Mood Tube is purposely designed for a softer, calmer lighting experience - ideal for use in maternity and neonatal environments.

Standard Sensory Lights vs Croyde Mood Tubes

Standard Sensory Tubes

- ✗ Bright, neon lights
- ✗ Cold, clinical colours
- ✗ Fast, flickering transitions

Croyde Mood Tubes

- ✓ 40% less light output
- ✓ Warm, calming colours
- ✓ Smooth, slow transitions

© 2020 CROYDE

CROYDE BIRTHING COUCH

Product Code: CBC



The innovative Croyde Birthing Couch is designed to support active and natural childbirth, helping women to find their most comfortable positions during labour and delivery. The Croyde Birthing Couch can be used as a full-length couch, folded over and used to support standing labour or used as a seat; and even be used as a double bed for post-natal resting.

Double covered in high quality, medical grade, breathable fabric, it is resistant to fluid ingress, easy to clean and compatible with hospital wipes and chlorine based cleaning solutions.

Mobile version coming soon. Contact us for more information.



CROYDE BIRTHING BEANBAG

Product Code: CBB

A barrel shaped beanbag designed to retain its shape and support upright positions. The bag can be shaped and moulded to the specific needs of the woman and is comfortable to cuddle and rest upon. Compatible with hospital wipes and chlorine based cleaning solutions.



CROYDE FLOOR MAT

Product Code: CFM

A large, comfortable floor mat delivering space for active movement during labour. It can be folded over to be stored in a compact, half-length position and is easily carried using the integrated handles. Compatible with hospital wipes and chlorine based cleaning solutions.



CROYDE POSITIONING PILLOW

Product Code: CPP

A triangular shaped beanbag that comfortably supports women when sitting upright or resting. Delivering a 60° backrest in the upright position and 30° in the lower position, it shapes and moulds to the body to enhance comfort whilst avoiding having to find lots of regular pillows. Compatible with hospital wipes and chlorine based cleaning solutions.



MULTITRAC BIRTH SYSTEM

Product Code: SDIDY



The Multitrac Birth System promotes close and active labours between the birthing person and partner. It is supplied with a plug-in Dullstein Birth Stool, Pelviball with fixed base, foot mat for comfortable kneeling and integrated rope system for vertical hanging. Multitrac is a complete system for supporting upright & active birth.

Rope Slings are supplied with the Multitrac as standard but are also available to purchase individually in a range of colours.



"Since purchasing the Multitrac earlier this year, women have been able to adopt positions that support physiological birth. The Multitrac enables the birth partner to provide support and comfort during labour and birth. It is easy to use, and the integrated rope sling has been a great tool for women."

Abigail Holmes, Deputy Head of Midwifery, Cardiff and Vale University Health Board

VARIOTRAC® SWIVEL SYSTEM

Product Code: VARIO D (Ceiling)/
VARIO S (Floor)/
VARIO W (Wall)



The Variotrac® swivel system is an integrated system that allows women to use an overhead rope in whichever position they like. It has 360° swivel capability and locks at 15° increments so you can position it above a birthing pool, couch, bed, birth stool, chair

or any other platform on demand. Rope slings are safely and easily attached onto any of the 5 hook positions and are available in a range of colours. Variotrac can be mounted into the ceiling, wall or floor to suit the room.

STATIC CEILING HOOK

Product Code: STAT-HAK

Our ceiling-mounted hooks can be installed into any birth room ceiling to give women the option to vertically hang, using a rope sling and adopt their most comfortable positions.



DULLSTEIN BIRTH STOOL

Product Code: DGH/ DGH-ALU



The Dullstein® Birth Stool is designed to support upright birth and is covered with a soft imitation leather cover to enhance comfort during usage. It has sealed seams to prevent fluid ingress and is available in a choice of colours – yellow, blue and orange with height options of 33cm or 38cm.

The Dullstein® can also be combined with the perfectly shaped inlay mat to support additional active birth positions. It is lightweight, easy to wipe clean and has a maximum safe working load of 135kg.



CROYDE KAYA BIRTH STOOL

Product Code: KBS

Strong, stable and easy to clean, this innovative birth stool is able to support a range of positions. It can be used under water for additional options in birthing pools and can also be used for rocking.



ERGOKNEELER

Product Code: ERGO001

The innovative Ergokneeler is an ergonomic kneeling system that protects knees and promotes good musculoskeletal health.

It can enhance comfort levels for Midwives who are supporting women during active labour and helps to reduce the compression forces and postural discomfort associated with low working height positions.



Moulded in polyurethane with active antibacterial properties, the ergonomic design offers supported kneeling and sitting positions whilst also being lightweight, portable and easy to clean.



CROYDE BIRTHING BALLS

Product Code: BBB(65/75)

Croyde birthing balls feature anti-burst technology and material that is 5 times thicker and more durable than standard gym balls.




PELVIBALL

Product Code: P-Ball (65/75)

The Pelviball can be used on its own or in conjunction with the Multitrac or Duoetto birthing systems.

The ball is attached to a base plate to offer maximum stability and security during labour, whilst also preventing it from inadvertent rolling. The ball is easy to wipe clean, is anti-burst and easy to top up with air over time.

Optional covers are available with disposable and washable versions.



HEIGHT ADJUSTABLE OVERBED COT

Product Code: BM1-HAC-75MM



The Croyde Height Adjustable Overbed Cot glides under the bed, putting the baby in front of Mum.

It allows constant and easy access, facilitating more frequent skin-to-skin contact, and promotes maternal bonding.

The mother can place the baby back in the cot without leaving the bed, reducing the risk of her falling asleep with the baby because the cot is too far away.

The Overbed Cot can be raised and lowered in one simple smooth movement to make access to the baby easier for women who have undergone caesarean sections or complicated births.

The open design makes it easy to clean, and a wide range of colours are available to suit any room theme.



Croyde Height Adjustable Cots are:



Used in over
80 NHS Trusts
across the UK

HEIGHT ADJUSTABLE STORAGE COT

Product Code: CAP-HACws-CSTD



Chosen by birth centres, postnatal areas, NICUs, SCBUs and HDUs, the Croyde Height Adjustable + Storage cot offers the convenience of customisable and practical under-cot storage.

The large cupboard is accessible from both sides of the cot with double soft-close doors, ensuring it is suitable for small rooms or tight spaces.

Height adjustment enables the cot to be raised or lowered to the desired height to limit bending or stooping, reducing musculoskeletal strain on the midwife's back and shoulders.

You can specify an optional foldable shelf or oxygen cylinder holder, and a wide range of colours are available to suit any room theme.



Croyde Height Adjustable Cots are:



Used in over
80 NHS Trusts
across the UK



AQUATHERM NEO WARMING MATTRESS

Product code: 7674-11282A73

The Croyde Aquatherm Warming Mattress is designed to provide controlled heat to newborn infants, particularly those who are premature or have difficulty maintaining their body temperature.

The mattress uses a fire retardant foam/ rubber and 'Aquagel' technology to deliver a non-water based warming system. A heating element is incorporated for safe and evenly spread heat as well as two highly-sensitive sensors for fail-safe temperature control.

When disconnected, the Aquagel technology retains residual heat to create a stable warming environment during transport.

The control unit has user-friendly settings to regulate the heat output of the mattress and also contains an alarm system for enhanced patient safety.



MEDWARM NEO WARMING MATTRESS

Product Code: IM-MS

The Medwarm Neonatal Warming Mattress sits in infant cots, offering controlled heat to keep newborns warm and comfortable - particularly premature babies and those that struggle to maintain their body temperature.

Medwarm stays fully operational when unplugged thanks to its rechargeable battery. The robust control unit will operate for up to four hours, keeping the mattress actively heated and reducing the risk of hypothermia.

The mattress operates at a temperature of 39°C, reaching 37°C in just 5-7 minutes, helping to protect an infant's core body temperature and avoid hypothermia after birth. The temperature is easily controlled, with an accuracy of 0.1°C.

The mattress core is made using memory foam, ensuring maximum comfort for baby whilst remaining lightweight for easy management and transport.



MAKE THIS A HOT COT

Add a Croyde Cot to create a fit for purpose Hot Cot. The Croyde Hot Cot provides a height adjustable, ergonomic and fit for purpose solution for neonatal warming.

Croyde Overbed Hot Cot with Aquatherm Neo - CAP-CHC01

Croyde Storage Hot Cot with Aquatherm Neo - CAP-CHCws-CSTD



MAKE THIS A HOT COT

Add a Croyde Cot to create a Hot Cot. The Croyde Hot Cot provides a height adjustable, ergonomic and fit for purpose solution for neonatal warming.

Croyde Overbed Hot Cot with Medwarm Neo - CAP-CHCM01 - M55

Croyde Storage Hot Cot with Medwarm Neo - CAP-CHCM03 - M55





HEIGHT ADJUSTABLE BABY ASSESSMENT TABLE

Standard Product Code: SPTM-02-D-HA
Compact Product Code: SPTM-02COMPACT-D-HA

Available in a standard or compact size, the Croyde Height Adjustable Baby Assessment Table allows clinicians to position the table at the perfect height for neonatal examinations and baby changing.

Depressing the foot pedal enables staff or parents to adjust the table to the ideal position and protect them from physical stress and strain. NIPE screening becomes instantly more comfortable and stress-free.

The high-height position makes cleaning easier, while the low-height position makes the table accessible for all heights of carers, including wheelchair users. It is lightweight and easy to move between bed spaces, with brakes to secure its location.

Choose from a storage cupboard with shelving or a drawer unit. Both units have soft-close doors and ample space for all items.



Standard	Compact
External Table Top Length (front to back) : 715mm	External Table Top Length (front to back) : 700mm
Internal Mattress Length (front to back): 650mm	Internal Mattress Length (front to back): 650mm
External Table Top Width (side to side): 810mm	External Table Top Width (side to side): 480mm
Internal Mattress Width (side to side): 680mm	Internal Mattress Width (side to side): 380mm
Minimum Height to Maximum Height: 810mm – 1100 mm	Minimum Height to Maximum Height: 810mm – 1100 mm



BABY CHANGE STATION

Product Code: CBS-M

The Croyde Baby Change Station has clinical application with a homely appearance. It's ideal for Birth Centres and Midwife- Led Units. The changing unit features a smooth tabletop surface with a curved upstand, deep, lockable drawers for safe storage, a soft oak finish and integrated 75mm castors with brakes for ease of movement.

The soft-close drawers provide enough space to store essential items. For example, you can store basic neonatal resuscitation equipment to replicate the home birth setup. Standardised equipment is then accessible, but not visible, a scenario that is well-received during diligence reviews.



CO-SLEEPER COT

Product Code: 150151-00 MED

The Croyde Co-Sleeper Cot is shaped to give immediate access to the baby with the safety of lipped edges and integrated side rails. It is designed for use in hospitals and helps to deliver a home-from-home environment.

CROYDE BEREAVEMENT SUITE BED

Product Code: XB1P8ABD106



CROYDE WALL BED

Product Codes:
Double wall bed: CDWB
Double wall bed & cabinetry: CDWC+C
Single wall bed: CSWB




The Croyde Electric Double Bed has been designed for hospital practicality with a non-clinical, homely appearance.

The head and foot boards are available in a range of colours to match your rooms, creating a non-clinical feel without compromising on practicality. The total weight capacity is 315kg and has integrated CPR release for emergency situations.

The backlit handset control allows the bed to be electrically adjusted for patient comfort, or trendelenburg position. The height of the bed can also be raised and lowered as needed.

The bed is supplied with a hospital specification mattress that is waterproof, fire retardant and breathable. It is also compliant with cleaning disinfectants, and fitted with a zip, protective flap and welded seam.



The Croyde Wall Bed is available as a single bed or a double bed and can be easily raised and lowered with just one hand.

The beds can be purchased independently or with alongside cabinetry such as wardrobes, cupboards, shelves or storage boxes.

Integrated lighting, USB charging points and power options are optional and a comprehensive delivery and installation service is available.

The cabinetry can be adapted to accommodate gas systems, nurse call systems or other fixed wall systems.

A wide range of colours and finishes is available and each bed is supplied with a hospital-specification mattress that fits perfectly to the frame and easily folds away with the bed.

All parts of the bed can be cleaned with hospital-approved wipes and chlorine-based solutions.



900MM
WIDTH

1200MM
WIDTH

1400MM
WIDTH



CROYDE BABY BATH

Product Code: CAP-TBB01

The Croyde Baby Bath is a height adjustable baby bath with out Easy-Pour Tub for quick and simple emptying.

The tub can be raised and lowered using the foot pedal to tailor the height to suit any care giver and improve ergonomics when using or emptying.

The Easy Pour Tub design features a pointed spout to make emptying the bath efficient whilst reducing spillages. The handle is easy to grip when tipping and can also be used as a small shelf to store soap sachets, cotton wool etc when in use.



TIPPER TUB

Product Code: TT-24000

This innovative mobile baby bath enhances caregiver ergonomics and eliminates the need to manually lift and carry heavy bowls of water.

The integrated lifting column allows the tub to be raised to the height of the sink or drain where it is then tilted without any

hand contact to empty without strain. The bowl is ergonomically shaped, heat-insulated to retain the warmth of the water and one-piece to make cleaning exceptionally easy.

Available in a choice of colours.

BREAST MILK FRIDGES

Product Code: BWR125

The Croyde Breast Milk Fridges are available in a range of sizes, from sixty litre capacity to four hundred litre capacity.

The digital control panel allows for the temperature selection, between 1-5 °C, with a single press 3 °C, setting for the optimum breast milk temperature. Audible and visual temperature alarms immediately alert healthcare practitioners of any major temperature fluctuations and help to ensure the quality of the milk.





CROYDE MATERNITY RECLINER

Product Code: ASTMR

A multi-use bedside package, perfect for maternity wards. The recliner facilitates ideal infant feeding positions for mother and baby. The reclining backrest allows babies to feed in the semi-recumbent chest-to-chest position that is most comfortable for the mother.

Padded armrests give comfortable arm placement and support the chosen feeding position. For partners and carers, the recliner positions them perfectly for skin-to-skin contact when bonding with the baby.

The padded foot section of the Croyde Maternity and Neonatal Recliner gently raises the level of the feet closer to the heart, improving lymphatic drainage, and reducing swelling in the feet, ankles, and legs after birth. By elevating the feet, the lower back also benefits from reduced strain.

The chair also allows overnight stays for partners without needing additional equipment, replicating the home environment and improving bonding with birth partners.



**UPRIGHT
BEDSIDE
CHAIR**

**RECLINE
FOR
FEEDING**

**RECLINE
FOR
SLEEPING**



CROYDE ELECTRIC MATERNITY RECLINER

Product Code: GeoSTIS20/GeoSTIS25

The Croyde Electric Maternity Chair aids sitting, reclining, and rising at the touch of a button, and supports women who have undergone invasive surgeries. The reclining backrest allows babies to feed in the semi-recumbent chest-to-chest position that is most comfortable for the mother. Padded armrests give comfortable arm placement and support the chosen feeding position.

It makes the immediate postnatal period more comfortable for women who have restricted mobility, who are experiencing pelvic girdle pain or are recovering from invasive procedures.

The Croyde Electric Maternity Recliner Chair facilitates overnight stays for partners without additional equipment needing to be brought into the ward. This helps to replicate the home environment and improve bonding with birth partners.



The Croyde Infant Feeding Pillow and Footstool are also available to purchase separately.

Available in a range of colours.



CROYDE KANGAROO RECLINER CHAIR

Product Code: TX15

The Kangaroo Recliner is an ideal choice for Maternity and Neonatal settings, combining full-recline comfort with smooth patient transfer functionality and a versatile selection of optional accessories to support both care and recovery.

The Croyde Kangaroo Recliner stands out with its unique ability to recline completely flat - a game-changer for exhausted parents in need of real, restorative rest.



Unlike many traditional recliners that can only partially recline - which can place unwanted pressure on sensitive areas, such as perineal or C-section wounds - this fully flat position provides true comfort and relief.

A fully reclined position is also ideal for caregivers to perform antenatal palpation and take symphysis fundal height (SFH) measurements.

Its ergonomic push handle allows caregivers to move patients smoothly and securely between clinical areas in an emergency, reducing strain and improving workflow.



CROYDE SLEEPER CHAIR

Product Code: STARGAZER

Perfect for use in hospital settings where overnight stays for partners/ dependants are required. The Croyde Sleeper Chair is easily deployed from chair to comfortable, full-length bed in just a few seconds. It also features a chaise longue position. Integrated wheels make moving simple and lightweight and it is designed to be wiped-down between uses in the hospital setting.

Available in a range of colours.



"The sleeper chairs are both spacious and comfortable, whilst not taking up too much room in our side rooms. They have helped us to ensure that our patient's relatives have a better experience during overnight stays, which are often in very difficult circumstances"

Inpatient Matron at Basildon & Thurrock, University Hospitals NHS Foundation Trust



Installed in over 70 NHS Trusts across the UK



MIMMAMÀ INFANT FEEDING CHAIR

Product Code: 365900 S/N

The MimmaMà is specifically designed for infant feeding to accommodate a variety of different positions including cradle, cross and rugby ball options.

The shape of the chair helps to maintain the correct position of the child on the nipple, avoiding fissures. It also allows the mother to alternate feeding positions, freeing the milk ducts and avoiding mastitis.

Ideal for twin feeding due to the dual backrest spaces. There is excellent support for the lower back, arms and shoulders.

The Croyde Infant Feeding Footrest and the Croyde Infant Feeding Pillow can be used alongside the chair to provide further support for mother and baby.

Available in a range of colours.



CROYDE FEATURE WALL PAPER

Product Code: Feature Wall (Wallpaper)

Bring colour to any hospital space using the Croyde Feature Hospital Wallpaper. A custom made wallpaper is an ideal way to cover a large area and make a great impact.

There is a wide range of designs available, from photographs of a local area to a calming coastal scene or landscape. These uplifting designs can transform clinical areas into welcoming spaces and help to create oxytocin-boosting birth environments.

Usually taking just a couple of hours to complete, this feature wallpaper can be installed without causing too much disruption to normal hospital routines and appointments.



CROYDE LIGHT CRYSTAL


Product Code: DOME

The Croyde Light Crystal cycles through a magical world of colours to deliver a calming, sensory experience.

It is placed in the corner of the room and patterns of light are projected through the crystal dome to produce a mesmerising, light display. The crystal works straight out of the box and just needs to be plugged in with no ongoing maintenance.

Perfect for areas striving to deliver relaxed environments.





QED BLANKET WARMING CABINETS

Product Code: BW-ECO

The QED Blanket Warming Cabinets are available in a range of sizes accommodating 4-50 cotton blankets in baskets. These cabinets can be used in range of settings, including birth environments, postnatal & neonatal, operating theatre recovery, critical care, A&E and other hospital wards where rapid patient warmth is important.

All QED Warming Cabinets are equipped with PID digital temperature controls. Blanket warmers have a maximum temperature of 80°C but can be set to a desired temperature.

The chamber is made from medical grade stainless steel and the external of the cabinet is powder painted to allow for easy cleaning.



QED FLUID WARMING CABINET

Product Code: FW-ECO

The QED Fluid Warmers are available in a range of sizes accommodating 9 to 240 one litre bottles of fluid.







Smaller cabinets are designed to fit underneath or on top of work surfaces, although stands are available if required. Larger cabinets are equipped with wheels for ease of movement, but we can also fit wheels on any of our cabinet sizes.

All QED Warming Cabinets are equipped with PID digital temperature controls.

Our fluid warmers arrive factory set to an optimal temperature of 37°C, but can be set to any desired temperature - based on the manufacturer's guidance of the fluid you're storing - up to a maximum of 60°C.



Indicative Sizing Options for Fluid and Blanket Warming Cabinets



MATERNITY CHAIR

Product Code: E600

Croyde maternity ward chairs are a fantastic addition to any unit and are ideal for use in post-natal wards and birth rooms alike. They have exceptionally comfortable visco-elastic foam seating to enhance the user's sitting/resting experience and a stylish, supportive design.

Available in a wide range of colours, they can be ordered to suit any room décor. The vinyl fabric is also infection control compliant and designed to be wiped down in-between use. Manufactured in the UK to exacting quality standards, the frames are stable, robust and provided with a 5 year warranty.



OVERBED TABLE

Product Code: F833

Croyde Overbed Tables are height adjustable and designed to allow the table to be used equally well with a bed or chair with its U-shaped base. They are lightweight, easy to move on wheels, durable for longevity of use and easy to clean. Table tops have lipped edges to prevent items sliding off and are available in blue, grey or beech finishes. Tilting or non-tilting top options are available.



VISITOR CHAIR

Product Code: PT10K/PT11K

Croyde visitor chairs are comfortable, practical and can be colour coordinated to patient chairs and lockers. They are available with a choice of fabric colours to complement modern room designs and colour schemes; whilst also having comfortable foam seating to enhance the user experience. Manufactured in the UK using only the highest quality materials, they are strong, durable and designed for use in the hospital setting.



CROYDE EXAMINATION COUCHES

Product Code: SM2560/ SM3560

A range of electric and hydraulically adjustable examination couches. All models have an excellent low height position to facilitate easy entry & exit, whilst also retaining a high height position that allows it to be set-up for taller clinicians and easy access for cleaning. Choose from 2 section or 3 section models and are available to order with retractable wheel bases or with castors is movement in-between rooms is required.

Available in a range of colours.



CROYDE SUTURING COUCH

Product Code: NV8683

The Croyde Suturing Couch features electric positioning for height, backrest angle and tilt position. It is supplied as standard with adjustable lithotomy leg rests, matching head support cushion, side support loops, debris tray and paper roll holder. An optional leg extension plate is also available. The 70cm platform width enhances comfort during usage and the large castors and centrally operated brakes make it lightweight and easy to transfer between rooms.



GRACIELLA GYNAE EXAMINATION CHAIR

Product Code: GRAC-CSS



The Graciella® has been designed with attention to detail, taking into consideration the time-consuming work of gynaecologists. With a unique two column design, the Graciella® enables easy positioning and increases patient trust before and during the examination.

COMFORTABLE POSITIONING

The positioning range allows the leg supports to be adjusted conveniently using the foot controller while ensuring that the patient remains comfortable.

There is also an integrated handset control which can be accessed from various angles due to the rotating mechanism.

FOOT PLATES

The unique foot plates enable comfortable and dignified leg positioning. Optional goepel leg rests are available and are easily interchanged.



PRE-SET POSITIONS

The Graciella can be easily adjusted to achieve pre-defined examination and entry/exit positions. Pressing one button co-ordinates multiple movements for examination and reverses them for exit afterwards.

CROYDE OPERATOR AND SADDLE STOOLS

Product Code: MC6164 / MC6144



Height adjustable operator stools and saddle stools. Comfortable seating with options for wheels, feet and brake-on seating.

Available in a range of colours.

CROYDE DIASYS EXAMINATION LIGHT

Product Code: DIASYS-MOB-S0-UK

A superb quality LED examination light capable of delivering high-level illumination for minor procedures.



Hands free brightness adjustment and automatic switch-on when you pull the light head into position, it is easy and practical to use when suturing.

Turns off after 30 minutes inactivity to avoid unnecessary usage.



CROYDE MINOR OPS E1 LIGHT

Product Code: MAR-LED E1

This examination lamp features high-quality, energy-efficient LEDs and a flexible extension arm to deliver illumination for minor procedures.

The mobile base of the light is stable, well balanced and durable with a radius footprint of just 350mm to make it easy to store.

The light head is available with a mobile base or can be mounted to a wall or ceiling.



All your maternity equipment needs, expertly supported.

Antenatal Equipment



Couches & Chairs



Examination Lights



Stools

Birth Room Equipment



Birthing Beds



Birthing Aids



Water Birth

Postnatal Equipment



Neonatal Warming



Cots & Baby Platforms



Phototherapy

Bereavement Suite Equipment



Double Beds



LED Sky Panels



Recliner Chairs



CROYDE MEDICAL



Scan to discover
the range