

COMBIPOOL BIRTH POOL

The Combipool is a beautifully shaped luxury birth pool that delivers an exceptional water birth experience.



The Combipool Birth Pool is supplied to an excellent specification, with an ergonomic built-in seat as well as integrated handles and foot supports for comfort and support.

Underwater LED lighting and an electronic filling system are also included.

The internal floor space is large enough to support freedom of movement whilst also offering a number of pre-defined positions that are easily achieved.

The shape and aesthetic of the pool is soft and gentle, contributing to a calming and inviting birthing environment.



BIRTHING SEAT & SUPPORTS

The integrated seat can be used to facilitate birth partner support whilst allowing the woman to be comfortable.

The shape of the pool allows the birthing person to sit comfortably with a reclined backrest and the perfectly positioned arm rests and foot plates help the woman to feel safe and secure.

The innovative grab handles can be used as hand grips or fixing points.



ELECTRONIC FILLING SYSTEM

The Combipool can be filled using a hospital supplied filling system or with its own, factory installed electronic filling system.

This system simply requires the pool to be connected to the utilities & saves valuable installation time and design cost.



UNDERWATER LED LIGHTS

Fully integrated under-water LED Lights illuminate the water for assessment/examination purposes, eliminating the need to submerge torches and compromise infection control.

These lights are also colour changeable so you can change the hue of the water to create an amazing sensory experience.



OPTIMUM HEIGHT & SHAPE

The sides are designed to achieve the optimum working height for the attending professionals. This minimises the risk associated with manual handling.

The shaped end of the pool also allows the midwife or birth partner to get closer to the water by placing their legs underneath the rim of the pool.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)



ENTRANCE STEPS WITH HANDRAIL

Using a unique 'flange' system, the optional steps securely lock into position for maximum stability.

This innovative design ensures that women are supported at every point during pool entry and exist as well as helping to avoid excessive twisting and uncomfortable leg lifting



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

