

CROYDE CLINNOVA SUTURING COUCH

Make suturing procedures easier with the Croyde Clinnova Suturing Couch.



The Croyde Clinnova Suturing Couch gives clinicians substantial access and flexibility. It helps women feel secure during procedures and examinations.

- Features include adjustable leg rests
- Adjustable leg rests
 - Crescent seat that folds backwards
 - Digital handset control
 - Removable debris tray, side support loops and paper roll holder
 - Optional foot section
 - Anti-bacterial powder coated frame

The Croyde Clinnova Suturing Couch is manufactured in the UK to a high specification.



ADJUSTABLE LEG RESTS

The lithotomy leg rests can be removed & reattached, and their height, rotation and angle easily adjusted. The fold-down armrests offer grabbing points to help women get into a comfortable position.



CRESCENT SHAPED SEAT

The crescent shape of the seat maximises perineal access. Below is a large stainless steel debris tray that captures fluid loss. It slides in and out and is easy to remove for cleaning.



DIGITAL HANDSET

The modern handset gives you digital control of the couch's position. The screen displays angles of the adjustable sections, such as the Trendelenburg tilt. When not in use, the plastic hook on the handset securely attaches it onto holders around the couch.



STORAGE & TRANSPORT

The couch has a footprint of 1.2m x 0.7m (without the optional foot section). It's ideal for storing in small areas when not in use. One person can move the couch between rooms by themselves due to the large wheels, lightweight frame and central braking system.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Choices](#)

[How to Buy](#)



DETACHABLE FOOT SECTION

We offer an optional detachable foot section that transforms the chair into a full-length examination platform. It's straightforward to remove and reattach whenever you need perineal access.



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

