

CROYDE HEIGHT ADJUSTABLE COTS

Chosen by birth centres, postnatal areas, NICUs, SCBUs and Neonatal Units, our height adjustable cots are available as an overbed option or with additional storage.



The Croyde Height Adjustable Cot range provides an ergonomic cot solution.

The height adjustability means that staff can position the cot at a safe working height to prevent musculoskeletal strain when performing bedside observations.

This also creates access for women with restricted movement or women who have undergone caesarean sections

A 6° head elevation tilt option supports neonatal digestive tracts in keeping milk down and moving wind up.



OVER BED OPTION

The Overbed Cot glides under the bed, putting baby in front of Mum, allowing constant and easy access, more frequent skin-to skin contact, and promotes maternal bonding.

The mother can place the baby back in the cot without having to leave the bed, lowering the risk of falling asleep with the baby because the cot is too far away.



STORAGE OPTIONS

Double door access from both sides mean that the cot can be used in small rooms or tight spaces.

The large cupboard provides storage for bulkier items and personal effects, to suit wards with limited storage space. Customisable shelves allow care equipment to be stored to suit the needs of the ward.



EASY TO CLEAN & MOVE

With its height adjustability, column design and open space construction, all areas are easily accessible for cleaning. It is compatible with hospital approved cleaning solutions and wipes and designed for use in high traffic wards. It is also lightweight and easy to move whilst still being stable and durable.



ERGONOMICS & ADJUSTABILITY

The cot can be raised and lowered in one simple movement to achieve the desired height to reduce musculoskeletal strain for the midwife. This also helps deliver effective new-born examination and enhanced observations at the bedside by adjusting to the height of the clinician and reducing maternal separation.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

CROYDE HEIGHT ADJUSTABLE OVERBED COT

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Options](#)

CROYDE HEIGHT ADJUSTABLE COT WITH STORAGE

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Options](#)



CROYDE MEDICAL
CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

