

# CROYDE POSITIONING PILLOW

Provide comfort for mothers in upright positions with the Croyde Positioning Pillow.



The Croyde® Positioning Pillow is a triangular beanbag with a supportive backrest.

It assists women in sitting upright by propping up their back, shoulders and neck. The inner beads help them feel supported in different positions during and post childbirth.

The pillow can support extra positions during labour and delivery. It also helps women who are experimenting with breastfeeding techniques and postnatal positioning.



## 60° BACKREST POSITION

When upright, the pillow's backrest elevation level is 60°.

The beads shape and mould around the back, neck, shoulders and head. This flexible movement provides optimum comfort.



## 30° BACKREST POSITION

In its secondary position, the pillow's backrest elevation level is 30°.

This shallower angle is suitable for resting or examination/palpation. The pillow moulds to individual body shapes and provides maximum comfort.



## CLEANING & INFECTION CONTROL

The Croyde Positioning Pillow is perfect for hospitals. It is made with high quality, medical-grade fabric that is durable, easy to clean & infection-control compliant.

You can wipe the pillow between uses with hospital wipes and chlorine-based solutions. The fabric is strong and resists fluid ingress. A double waterfall protection flap covers the zip and PU tape seals and protects all seams.



## ACTIVE BIRTH PRODUCTS

The Croyde Positioning Pillow perfectly complements other products in the Croyde Active Birth collection, helping midwives to maximise birth position possibilities.

The Croyde Birthing Beanbag, Croyde Floor Mat and Croyde Birthing Couch can be used alongside the Croyde Positioning Pillow to create a range of positioning options.

## EVERYTHING YOU NEED TO KNOW

Click the links below to access the

Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Videos](#)



## BIRTHING BED COMPATABILITY

You can use the pillow with a delivery bed to support advanced and active positioning,

The inner beads also enhance comfort and can provide support to help make long labours more manageable.



CROYDE MEDICAL

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[www.croydemedical.co.uk](http://www.croydemedical.co.uk)

Croyde Medical, Unit A4 GBP1, Grove Business Park,  
Downsview Road, Wantage, OX12 9FA

