

CROYDE TIPPER TUB

The Croyde Tipper Tub is an innovative ergonomic mobile baby bath.



The Croyde Tipper Tub enhances caregiver ergonomics by eliminating the need to manually lift and carry heavy bowls of water.

The lifting column allows the tub to be raised to the height of the sink or drain where it is then tilted and emptied without strain.

The tub is easy to clean with a one-piece finish free of any joints or edges. It is also available in a range of colours.



ERGONOMICS

The Croyde Tipper Tub central column can be adjusted to fit the height of a sink where it would need to be emptied.

The head can be tipped to allow water to exit the tub without the need of tipping the entire bath.

These design choices should reduce strain on staff when bathing newborns.



COLOUR RANGE

The Croyde Tipper Tub is available in a wide range of colours to suit any hospital ward.

The most popular colour is Signal White RAL 9003 but there are over 40 other colours to choose from.



CLEANING

The Croyde Tipper Tub is exceptionally easy to clean with a one-piece finish that is free of any joint or edges. It makes wiping and drying the tub a simple task.

The adjustable lifting column and tilt handle are made of stainless steel for ease of cleaning.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the
Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Options](#)

[CLICK HERE FOR INFORMATION ON HOW
TO BUY THE CROYDE TIPPER TUB](#)



CROYDE MEDICAL

Maternity Equipment Specialists

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

[**www.croydemedical.co.uk**](http://www.croydemedical.co.uk)

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX

