

GENEROSA BIRTH POOL

The Combipool is a beautifully shaped luxury birth pool that delivers an exceptional water birth experience.



Designed with an ergonomic built-in seat as well as integrated handles and foot supports for comfort and support, the Generosa Birth Pool creates an optimal environment for the birthing person.

The compact footprint gives the same benefits as other birthing pools, whilst making it accessible for areas where floor space is limited.

The standard specification birth pool consists of a white inner and outer shell, integrated birth seat and arm rests and an integrated pop-up waste with external release.



BIRTHING SEAT & FREEDOM OF MOVEMENT

As standard, the Generosa is supplied to an excellent specification, featuring an integral birthing seat and handles.

This pool allows the birthing person to move freely without restriction before and during birth, allowing them to find their most comfortable position.



INTEGRATED POP-UP WASTE

The integrated pop-up waste is controlled by simply rotating the dial to open the drain.

The dial is conveniently & ergonomically located in the external shell of the pool to make emptying the pool an easily accessible & simple procedure.



ENTRANCE STEPS WITH HANDRAIL

Using a unique 'flange' system, the steps securely lock into position for maximum stability.

This innovative design ensures that women are supported at every point during pool entry and exist as well as helping to avoid excessive twisting and uncomfortable leg lifting



UNDERWATER LED LIGHTS

Fully integrated under-water LED Lights illuminate the water for assessment/ examination purposes, eliminating the need to submerge torches and compromise infection control.

The variety of coloured light options provide mood lighting.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)



INTEGRATED FILLING SYSTEM

The optional, fully integrated filling system gives nurses and midwives control over each feature of the birth pool as well as providing access for servicing.

This includes a thermostatic mixing valve to control water temperature, autofill button, light controls and shower controls.



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

