

LIFETHERM 2000 WARMING BED

From first breath to full stability - Lifetherm 2000 ensures every newborn starts life warm, safe and strong.



The Lifetherm 2003 Warming Bed combines neonatal warming with critical care support.

It ensures safe, reliable thermoregulation while streamlining workflows and improving outcomes for vulnerable infants.

With both underbody warming (precisely controlled between 28–38.5°C) and overhead radiant warmth, heat support can be gradually reduced, allowing babies to develop their own thermoregulation safely and comfortably.



REDUCED INCUBATOR SEPARATION

It supports early stabilisation by delivering precise thermoregulation, reducing time spent separated from parents in an incubator.

The integrated overhead radiant warmer provides immediate thermal support during critical interventions. An optional lid is available to further improve heat retention.



STREAMLINING CLINICAL ACCESS

The Lifetherm 2003 reduces care pressures by ensuring effortless access and practical support at all times.

Hinged side panels allow immediate access to the baby for both parents and staff, while sturdy rails accommodate essential accessories such as trays, ventilators and Apgar timers.



OPTIMISED MANUAL HANDLING

Safe, efficient handling is essential for both staff and infants - and a key focus of the CQC in maintaining care standards. Quiet, vibration-free height adjustment enables ergonomic and strain-free positioning for staff, without causing unnecessary disruption to the infant.



VERSATILE STORAGE SOLUTIONS

An optional storage drawer block can be configured on the front or back of the lifting column, offering convenient storage access to staff. This block features two drawers and a shelf, all of which can be independently rotated to provide easy, versatile access.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the

Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

