

# MEDWARM NEONATAL WARMING MATTRESS


The Medwarm Neonatal Warming Mattress sits in infant cots, offering controlled heat to keep newborns warm and comfortable.



The Medwarm Neonatal Mattress provides under-body heat for newborn infants to stay warm and maintain stable temperatures in their cot during transfers or procedures, helping to prevent hypothermia and related complications.

The mattress operates at a temperature of 30-39°C. Healthcare providers can easily control the temperature of the mattress with an accuracy of 0.1°C.

The mattress core is made using memory foam, ensuring maximum comfort for baby whilst remaining lightweight for easy management and transport.



## SUSTAINABLE BY NATURE

Designed to be reusable, the Medwarm Neonatal Mattress offers sustainable warming. The mattress is covered with Polyurethane - waterproofing it against all liquids and making the surface easy to clean and disinfect before its next use.



## SAFETY PRIORITISED

Built-in alarms will detect and alert if the mattress gets too hot or in case of malfunction, disconnection or short-circuit. A carbon fibre layer allows X-rays to take place without moving the baby, ensuring continuous warmth and comfort during imaging.



## ALWAYS READY

The mattress pre-warms and stays at the set temperature without needing weight, so babies benefit from immediate warmth when placed on it. The mattress will then maintain its warmth the whole time it is turned on, ensuring consistent warming when its needed.



## BATTERY OPERATION

Unlike other heated mattresses, the Medwarm stays fully operational when unplugged thanks to its rechargeable battery. The robust control unit will operate for up to three hours, keeping the mattress actively heated and reducing the risk of hypothermia - even in cold or drafty environments.

This facilitates transfers, offers a theatre-ready design and allows for enhanced parental experiences.

## EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Warranty & Servicing](#)

[Cleaning Instructions](#)



CROYDE MEDICAL

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[www.croydemedical.co.uk](http://www.croydemedical.co.uk)

Croyde Medical, Unit A4 GBP1, Grove Business Park,  
Downsview Road, Wantage, OX12 9FA