

# MEDWARM PATIENT WARMING SYSTEM

Offering controlled heat to help patients maintain comfortable body temperatures in theatre rooms, reducing the risk of hypothermia.



The Medwarm Patient Warming System is made up of the Medwarm Mattresses and Control Unit as standard, although can also be used with the Medwarm Blankets.

It offers efficient and intuitive conductive heating, ensuring patients' body temperatures are effectively regulated.

The mattress operates a temperature range of 30-40°C, which is adjustable in 0.1°C increments for precise control and accuracy. Heat is then distributed evenly and consistently throughout the mattress.



## SIMPLE & DURABLE DESIGN

The Medwarm W-300 Control Unit offers a perfect blend of simplicity and durability. With a user-friendly interface, the W-300 features a modern, robust design built to withstand frequent use and the demanding conditions of hospitals.



## SUSTAINABLE SOLUTION

The Medwarm Mattress provides a reusable and eco-friendly solution for patient warming that doesn't use consumables.

With a 10-year lifespan, it offers extensive usage before replacement is required - offering unmatched longevity.



## INFECTION CONTROL

The Medwarm Mattress is perfect for environments where infection control is critical.

It uses conductive heat - eliminating air movement and minimising the risk of disturbing sterile fields or spreading airborne contaminants.

## COMFORT & SAFETY PRIORITISED

Built-in alarms make up a comprehensive security system that detects and alerts if the mattress gets too hot, malfunctions, disconnects or short-circuits.

The mattress uses X-ray-permeable carbon fibre, allowing patient warming during scans, while memory foam adds comfort and helps prevent bed sores during long procedures.



## EVERYTHING YOU NEED TO KNOW

Click the links below to access  
the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Warranty & Servicing](#)

[Cleaning Instructions](#)



**CROYDE MEDICAL**  
Maternity Equipment Specialists

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[\*\*www.croydemedical.co.uk\*\*](http://www.croydemedical.co.uk)

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX