



CROYDE MEDICAL

MULTITRAC BIRTHING SYSTEM

Support partner engagement and active birth with the Multitrac Birthing System.



The Multitrac Birthing System promotes close and active birth between mother and partner.

The seat is fitted with a comfortable padded seat and backrest which is available in three colour options of yellow, apricot or blue.

The Multitrac is also easily compatible with plug-in birth aids to provide various positioning options.

It also comes with an additional integrated rope sling – useful for birthing rooms without permanent, ceiling-mounted rope slings.



BIRTH PARTNER ENGAGEMENT

The Multitrac encourages a close and active birth. As the woman uses a plug-in birth aid, her partner can sit on the chair behind, and be at the same height. The woman can also turn around and kneel on the mat whilst leaning onto a birth aid or her partner's lap.



PLUG-IN BIRTH AIDS

The Multitrac is extremely flexible and modular in its design. The base plate is compatible with the Pelviball and the Dullstein birth stools. Both can easily plug into the system or be used on their own.



FLOOR MATS

Optional inlay mats and floor mats further promote active birth. These deliver a more comfortable surface when kneeling around the seat, Pelviball or Dullstein birth stool.

The mats are waterproof and fitted with welded seams. They are easy to wipe down between uses.



MOVING AND HANDLING

The Multitrac Birth System is sturdy, supportive and highly stable during usage.

One person can simply move the system by tipping it back 45 degrees. Integrated castors on the rear wheels make it easy to move around and clean.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Cplour Options](#)

[CLICK HERE FOR INFORMATION ON HOW TO BUY](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

