



CROYDE MEDICAL

NESTIE MATERNITY BED

The Nestie Maternity Bed is suitable for use during intrapartum and postnatal care, as well as for facilitating partner stays.



The Croyde Nestie Maternity Bed has been designed with both maternity staff and the birthing person in mind, to create a home-from-home experience in a clinical setting.

The Nestie is available with a choice of different colour head and footboards, is available in single and double-width options, and also comes with detachable lithotomy poles to assist with advanced positioning during labour or for suturing.



ELECTRIC ADJUSTABILITY & ONBOARD CPR

The built-in handset can be used to control the height of the bed and adjust the backrest, thigh section and Trendelenburg tilt.

The auto-chair function allows the bed to be moved into a chair position without having to manually adjust individual settings.

The integrated CPR ensures the bed is ready to accommodate emergency scenarios when required.

WIDTH & WEIGHT CAPACITY

The Nestie Maternity Bed is available as a single bed or in one of three double sizes; 1070mm, 1200mm or 1370mm wide. This gives options to allow for more space to an individual occupant or enough space for two people to stay together comfortably.

The bed has a safe working load of 400kg which may also be of interest for the care of bariatric maternity patients.



LITHOTOMY POSITIONING

The bottom of the bed is equipped with a medirail that allows lithotomy poles to be attached to the bed. The poles can be used to support lithotomy positioning or foot plate positioning to support advanced birth positions or suturing as required.

The positioning pillow allows the woman to be comfortably placed in a semi-recumbent position during these procedures.

The under bed clearance comfortably permits the suturing stool and midwife to get closer to the perineum during suturing.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the

Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[CLICK HERE FOR INFORMATION ON HOW TO BUY](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

