

OCEANIS BIRTH POOL

The Oceanis Birth Pool is a warm, spacious birth pool delivering freedom of movement in the water.

Manufactured to the highest quality standards, the Oceanis pool looks inviting and feels luxurious.

Designed to support women and midwives to achieve their natural childbirth objectives, the Oceanis has a wealth of unique features and customisable design options to complement any delivery room.

The standard specification birth pool consists of a white inner and outer shell, integrated birth seat and arm rests and an integrated pop-up waste with external release.



BIRTHING SEAT & FREEDOM OF MOVEMENT

The large, spacious base of the pool allows women to move freely without restriction and the integral birthing seat & perineal cut-out with armrests create a comfortable environment.

The armrests and footrest positioning plates also offer additional options for women to find their most comfortable position.



INTEGRATED POP-UP WASTE

The integrated pop-up waste is controlled by simply rotating the dial to open the drain.

The dial is conveniently & ergonomically located in the external shell of the pool to make emptying the pool an easily accessible & simple procedure.



ENTRANCE STEPS WITH HANDRAIL

Using a unique 'flange' system, the steps securely lock into position for maximum stability.

This innovative design ensures that women are supported at every point during pool entry and exist as well as helping to avoid excessive twisting and uncomfortable leg lifting



UNDERWATER LED LIGHTS

Fully integrated under-water LED Lights illuminate the water for assessment/ examination purposes, eliminating the need to submerge torches and compromise infection control.

The variety of coloured light options provide mood lighting.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)



INTEGRATED FILLING SYSTEM

The optional, fully integrated filling system gives nurses and midwives control over each feature of the birth pool as well as providing access for servicing.

This includes a thermostatic mixing valve to control water temperature, autofill button, light controls and shower controls.



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

