

# QED COMBI WARMING CABINETS

A selection of combined blanket/fluid warming cabinets in a range of sizes to suit your needs..



QED Combi Warming Cabinets can be used in various areas, including:

- Birth Environments
- Postnatal Care
- Neonatal Care
- Operating Theatre Recovery
- Critical Care
- A&E
- And other hospital wards where warm blankets and fluids need to be co-located.

Optimised for warmth, they feature a quiet internal fan to ensure even, warm air distribution around the chamber.



**5 YEAR WARRANTY**

## LEADING WARRANTY

All QED Warming Cabinets are supplied with an **industry-leading five-year warranty**.

Labour is included for the first year, while spare parts are supplied for the **full five-year period** to reduce maintenance costs.



## TEMPERATURE CONTROL

Blanket warmers have a maximum temperature of 80°C and fluid warmers have a maximum temperature of 60°C.

Both can be set to a desired temperature using the PID digital temperature controls. Overheat shutoff is also included.



## EASY CLEANING

Each cabinet is fitted with a stainless steel interior for easy cleaning.

The stainless steel interior and powder-coated exterior are compatible with hospital-approved wipes and chlorine-based cleaning solutions.



## VAST SIZE RANGE

Combi Warmers come in a wide range of sizes dictated by the individual size of the fluid warmer and the blanket warmer.

These sizes do not have to be equal, meaning your chosen Combi Warmer can have more space dedicated to the fluid warmer or blanket warmer.

Fluid warmers compatible with the combi warmers range hold 15 to 72 one litre bottles of fluid.

Blanket warmers compatible with the combi warmers range hold 4 to 24 blankets.

## EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Warranty & Servicing](#)

[Cleaning Instructions](#)



**CROYDE MEDICAL**

Maternity Equipment Specialists

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[\*\*www.croydemedical.co.uk\*\*](http://www.croydemedical.co.uk)

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX