

QED WALL-MOUNTED CRT CABINETS

The QED Wall-Mounted Medicine CRT Cabinets act as traditional wall-mounted medicine cabinets whilst offering controlled room temperature conditions.



The QED Wall-Mounted Medicine CRT Cabinets are designed around traditional wall mounted medicine cabinets, yet as a space-saving alternative to free-standing CRT units - all whilst offering a higher level of security.

They feature integrated temperature control, ensuring medications are stored in room temperature conditions to maintain their efficacy and lifespan.



REDUCING MEDICINAL WASTE

The cabinets feature integrated temperature control to keep medicines below 25°C, which is crucial for maintaining efficacy in warmer weather.

By ensuring proper storage, the QED MED-CRT cabinets help prevent medicine degradation and significantly reduce waste.



LEGISLATION COMPLIANT

Fully compliant with all current legislation for medicine storage in healthcare settings, each cabinet is constructed from 1.2mm thick steel with a fully welded design for durability.

Anti-jimmy folded returns and a lock resistant to picking and drilling are also included.



TEMPERATURE CONTROL

Each unit features a pharmacy-grade digital temperature controller that is set to control between 15°C and 20°C.

There are also built-in alarms for added safety and monitoring should temperatures fall outside of this range.



CD CABINET COMPATIBLE

If needed, the QED MED-CRT cabinets can be fitted with a secure internal controlled drugs (CD) unit.

This provides a convenient and compliant solution for storing controlled substances safely within the main unit whilst keeping essential medications organised and easily accessible.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Warranty & Servicing](#)

[Cleaning Instructions](#)



CROYDE MEDICAL

Maternity Equipment Specialists

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

[**www.croydemedical.co.uk**](http://www.croydemedical.co.uk)

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX