



CROYDE MEDICAL  
Maternity Equipment Specialists

# QED PHARMACY REFRIGERATORS

QED Pharmacy Fridges give you reliable, secure storage for critical medicines and vaccines - so you can focus on patient care with total confidence.



Designed with convenience and reliability in mind, the QED Pharmacy Fridges feature automatic defrost and condensate evaporation systems that enhance both performance and ease of use. Together, these features help maintain an optimal storage environment while saving time and effort for busy healthcare professionals.

They are engineered for precision and reliability, being factory-set to a stable 5 °C within a controlled range of 2–8 °C. This ensures all medicines are kept at an optimal temperature - maintaining their efficacy and shelf life.



## UNRIVALLED WARRANTY

With an industry-leading **five-year warranty**, the QED Pharmacy Refrigerators offer unmatched peace of mind and long-term reliability.

They are fully CE-certified and compliant with all current refrigeration guidelines.



## PRECISION CONTROL

Advanced digital microprocessors ensure accurate temperature regulation, while built-in internal fans provide rapid cooling recovery after every door opening.

The fridges also automatically log min. and max. temperatures, and feature high/low alarm alerts.



## ROBUST SECURITY

Lockable doors provide robust security for safeguarding high-value vaccines and controlled medications around the clock.

There are also door-open alerts and power-failure notifications to ensure critical stock is always protected.

PR300



PR30



PR60



PR500



PR640



PR900



## EXTENSIVE MODEL RANGE

The QED Pharmacy Refrigerators are available in 11 versatile sizes, from space-saving under-counter models to high-capacity freestanding units ranging from 30 L to 1,390 L.

With options for solid or glazed doors and convenient bench-top designs, there's a solution for every setting - whether you're outfitting a compact clinic, a busy hospital or a clinical trial facility demanding reliable, scalable cold storage.

## EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Warranty & Servicing](#)

[Cleaning Instructions](#)



**CROYDE MEDICAL**

Maternity Equipment Specialists

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[www.croydemedical.co.uk](http://www.croydemedical.co.uk)

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX