


CROYDE WALLBED

A beautiful fold-down bed for use in hospitals



The Croyde Wall Bed is available as a single bed, or a double bed to suit your needs. The bed is easy for one person to fold down and put away.

Each bed is also available with a range of alongside cabinetry such as cupboards, drawers, wardrobes or top-boxes to accommodate storage requirements.

The bed is supplied with a hospital-specification mattress that can be cleaned using hospital approved cleaning methods and is designed for repeated usage in high-throughput wards.



BESPOKE CONFIGURATIONS

The Croyde Wall Bed is supplied as a single bed or a double bed in its standard configuration. Bespoke designs can be created to fit your space and requirements. Cupboards, wardrobes and top-boxes can be made to fit a specific space with floor-to-ceiling coverage. Cabinets can also be adapted to accommodate gas supplies or nurse call systems. Integrated power sockets and lights can also be included.



HOSPITAL SPECIFICATION MATTRESS

The infection control compliant mattress fits in the frame and folds away with the bed. It can be cleaned with hospital-approved disinfectants and wipes.

The mattress is also comfortable for sleeping, whilst also having a breathable, two-way stretch, moisture vapour permeable cover, with welded seams and a covered zip for internal access.



WALLBED - MIX

The Mix Wallbed has a spring assisted open and close mechanism for smooth and easy operation.



WALLBED - SPACE

The Space Wallbed has a piston based mechanism for an easy open and close

WALLBED - MIX EVERYTHING YOU NEED TO KNOW

Click the links below to access the
Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour & Design Options](#)

[How to Buy](#)

WALLBED - SPACE EVERYTHING YOU NEED TO KNOW

Click the links below to access the
Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour & Design Options](#)

[How to Buy](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

[**www.croydemedical.co.uk**](http://www.croydemedical.co.uk)

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

