



Croyde Variable Height Cots

Safer, closer Mum to Baby positioning



- Variable height neonatal cots with a clean, contemporary appearance
- Low height frame enables under-bed positioning of the cot
- Allows baby to be closer to the mother when in bed for skin to skin contact
- Available in **any** colour to match furniture / room décor
- Exceptionally easy to clean with excellent access to all areas
- Trendelenburg tilt +12° to permit head elevation
- Lightweight and easy to transport with 75mm castors
- Brake castors on wheels for added security / stability
- Colour co-ordinated bedside cabinets also available



CROYDE MEDICAL

Maternity Equipment Specialists

- Standard specification includes choice of **any colour**, storage basket on the column, manual Trendelenburg adjustment, hi-lo variable height adjustment, brakes with castors, clear plastic crib with mattress
- Clear plastic crib easily removed for cleaning
- ABS plastic cover protects base frame
- 1 year warranty
- Optional medi-rail can be mounted on the column of the cot to hang useful accessories such as warming mattress pumps
- Optional additional storage basket on the inside of the column



Specification

External cot length	780mm
External cot width	480mm
External cot height	970mm
Height adjustment range	970mm – 1240mm
Storage Basket Height	335mm
Storage basket width	405mm
Storage basket height	90mm
Internal crib length	660mm
Internal crib width	410mm
Internal crib height	220mm
Trendelenburg	+12°
Lead Times	4–6 weeks



CROYDE MEDICAL

Maternity Equipment Specialists

Croyde Medical, Unit 7, Regal Way,
Faringdon, Oxon SN7 7BX