

Diasys Procedural Light

High quality light for examination & suturing



- Motion-Matic sensor delivers intelligent illumination
- The Diasys automatically turns on when moved
- No on/off switch saves time and energy
- Automatically shuts off after 30 mins of in-activity
- Point the light head upwards to turn off
- Motion-Matic design eliminates the possibility of the light being left on and lamps burning-out
- Compact storage base & easy to transport
- Integrated push handle for ergonomic movement



CROYDE MEDICAL

Maternity Equipment Specialists

- Handsfree adjustment of brightness
- Simply wave the hand past the sensor to change the light intensity
- Ideal for when the clinician is gloved-up
- 4 * LED lights deliver excellent illumination
- Excellent range of adjustment
- Spring arm is stable & holds position
- Sterilizable handle
- Easy to wipe clean
- Optional battery back-up
- Excellent shock absorption
- No yellowing over time
- 1 year warranty

Available in mobile, wall or ceiling mounted configurations

Specification

Light source	LED
Maximum central lighting	65,000 LUX
Adjustable luminous intensity	30,000 to 65,000 LUX
Fixed field diameter	14cm
Depth of field	128cm
Colour rendering index	96
Red colouring rendering	95
Colour temperature	4200K
LED Life time	60,000H
Max consumption	7 W
Light base dimensions	30cm * 24.5cm



CROYDE MEDICAL

Maternity Equipment Specialists

Croyde Medical, Unit 7, Regal Way,
Faringdon, Oxon SN7 7BX