

Croyde Bedside Furniture

Maternity Ward High Back Chairs



SKU – E601

Croyde Maternity Ward High Back Chairs have a modern, fresh design to help enhance the look and feel of any maternity unit.

They are supplied with visco-elastic memory foam cushioning to deliver exceptional levels of comfort, housekeeping wheels for ease of movement within the room and ergonomic arm rests to support patient standing.

With strong, sturdy frames, 18" seat height position, wipe-clean vinyl upholstery and the choice of a wide range of colours, these patient chairs are perfect for maternity wards.



CROYDE MEDICAL

Maternity Equipment Specialists

Croyde Bedside Furniture

Maternity Ward Visitor Chairs



SKU – F833



SKU – F835



SKU – V020



SKU – V300



SKU – V310

Croyde maternity ward visitor chairs are supplied with comfortable foam seating to enhance the user experience and are available in a choice of colours to complement any room scheme.

Wipe-clean vinyl upholstery ensures ease of cleaning and infection control compliance; whilst all frames are strong, lightweight and stackable.

Supplied as standard with 18" seat height positions for optimal ergonomics and options for low back or medium back supports, Croyde visitor chairs are perfectly tailorable to any unit.



CROYDE MEDICAL

Maternity Equipment Specialists

Croyde Bedside Furniture

Single Piece Visitor Chairs



SKU -T110AN - GREY

Croyde single piece visitor chairs are strong, lightweight & durable. They are available in a choice of contemporary, vibrant colours and offer a fresh, modern appearance

The single moulded form makes it quick and simple for the chair to be wiped down and all options are stackable. Arm rests are available with all colour options.



SKU -T100AN- MOSS GREEN



SKU - T110 - YEL



SKU - T100AN-SKY BLUE



CROYDE MEDICAL

Maternity Equipment Specialists

Croyde Bedside Furniture

Bedside Lockers & Overbed Tables



SKU – VISA4+YEL



SKU – VISA4+Colour

Croyde bedside lockers offer significant bed space storage with two drawers (top drawer lockable) and one large cupboard. They are available in a range of colours and are also equipped with wheels for ease of movement. Designed to be cleaned in hospital wards.

Croyde overbed tables are gas spring height adjustable and ergonomically designed with a U shaped base to allow the table be used equally well with a bed or chair. Table tops have lipped edges to keep items safe on the table top and are available in a choice of blue, grey or beech finishes. These tables are lightweight, durable and easy to clean



SKU – PT10KCMBC



Croyde Bedside Furniture

Dimensions and Specifications

	Croyde Maternity High Back Chair E601	Croyde Visitor Chair F833 / F835	Croyde Visitor Chair V300	Croyde Visitor Chair V300
Seat height	482mm	483mm	457mm	457mm
Seat width	502mm	500mm	483mm	483mm
Seat depth	457mm	570mm	482mm	482mm
Overall height	1160mm	850mm	864mm	1041mm
Overall width	620mm	550mm	584mm	610mm
Overall depth	457mm	530mm	610mm	660mm
Safe Working Load	127kg	127kg	127kg	127kg
Crib	5	5	5	5
Warranty	5 years frame 1 year fabric	5 years frame 1 year fabric	5 years frame 1 year fabric	5 years frame 1 year fabric
	Croyde Visitor Chair T100AN	Croyde Visitor Chair T110AN	Croyde Overbed Table PT10KCMBC	Croyde Bedside Locker VISA4
Seat height	460mm	460mm	N/A	N/A
Seat width	480mm	480mm	N/A	N/A
Seat depth	465mm	500mm	N/A	N/A
Overall height	897mm	790mm	710mm > 1080mm	850mm
Overall width	635mm	500mm	900mm (Top)	462mm
Overall depth	534mm	514mm	420mm (Top)	422mm
Safe Working Load	158kg	158kg	25kg	N/A
Crib	N/A	N/A	N/A	N/A
Warranty	5 years	5 years	1 year	5 years

